NOTE: This page is not to be printed as part of the print job.

CKMB Calibrator

Contents

REF	Contents	
09318028	2 vials of lyophilized Low Calibrator	CAL L
(672174)	2 vials of lyophilized High Calibrator	CAL H
	Calibrator Assigned Value card	
	barcode labels	

10994373 Rev. B, 2014-09

Intended Use

For *in vitro* diagnostic use in calibrating ADVIA Centaur® systems CKMB assays.

Calibrator Description

Refer to the Calibrator Assigned Value Card for the assigned value of each calibrator. The assigned value is traceable to the standardization of the CKMB assay. For additional information, refer to the reagent instructions for use.

Volume	Ingredients	Storage	Stability
2.0 mL/vial	After reconstitution, low or high levels of CKMB in equine serum with gentamicin sulfate, preservatives, and protein stabilizers	2–8°C	Lyophilized—until the expiration date on the vial label or reconstituted—14 days or onboard—4 hours

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

For *in vitro* diagnostic use.

Preparing the Calibrators

1. Add 2.0 mL of reagent water that is 20–25°C into each calibrator vial using a volumetric or precision pipette.
- NOTE:** For information about reagent water, refer to the system operating instructions.
2. Gently mix the vials for 30 minutes on a mechanical mixing device (for example, a rocker or a rotator) at room temperature (20–30°C) to allow the lyophilized material to dissolve.
3. Gently swirl and invert the vials until homogeneous.

Using the Calibrator Assigned Value Card

Each lot of calibrators contains a Calibrator Assigned Value card to facilitate entering the calibrator values on the system. Enter the values using the barcode wand or the keyboard. For more information about entering calibrator values, refer to the system operating instructions or to the online help system.

Using the Barcode Labels

NOTE: Calibrator barcode labels are lot number specific. Do not use barcode labels from one lot of calibrators with any other lot of calibrators.

Affix the Low and High Calibrator barcode labels to the calibrator sample cups so that the system recognizes the sample as a calibrator. Place the barcode label on the sample cup with the readable characters oriented vertically.

Scheduling the Calibrators

Refer to *Calibration Interval* in the appropriate product insert for the frequency of calibration.

1. Schedule the calibrators to the workload:
 - a. Schedule CKMB Calibrator for the CKMB test.
 - b. Continue to schedule quality control and patient samples for the required tests.
2. Label the sample cups, one with a Low Calibrator barcode label and another with a High Calibrator barcode label.
Because calibrator barcode labels are lot number specific, use barcode labels that correspond to the lot number of the calibrator used.
3. Gently mix the Low and High Calibrators.
4. Dispense at least 1.0 mL of Low and High Calibrator into the labeled sample cups.
The sample cups are marked at 1.0 mL intervals to assist in filling.
The required calibrator volume depends on how many assays are being calibrated using these calibrators and the number of calibrator replicates.
5. Load the calibrator sample cups on the system.
The Low Calibrator cup position must precede the High Calibrator cup position.
Ensure that the CKMB reagents are loaded on the system.
6. Continue loading quality control and patient samples.
7. Start the system, if required.

Limitations

- Ensure that the water used to reconstitute CKMB Calibrators is at a temperature of 20–25°C. Using water with a temperature less than 20°C may result in higher recoveries. Using water with a temperature greater than 25°C may result in lower recoveries.
- Do not return any calibrators back into the vials after calibration because evaporation can occur, which may affect performance.
- Dispose of any calibrator remaining in the sample cups after 4 hours.
- Do not refill calibrator sample cups when the contents are depleted. If required, dispense fresh calibrators.

Disposal

Dispose of hazardous and biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with all federal, state, and local requirements.

Technical Assistance

For customer support, please contact your local technical support provider or distributor.

www.siemens.com/diagnostics

ADVIA Centaur is a trademark of Siemens Healthcare Diagnostics.

© 2014 Siemens Healthcare Diagnostics. All rights reserved.

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Siemens Healthcare Diagnostics Inc.
Siemens Healthcare Diagnostics Ltd.
Siemens AG
Diagnostics Inc.