Toxoplasma IgM (μ -Capture)

**For the Qualitative Detection of
IgM Antibodies to *Toxoplasma gondii***

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Toxoplasma IgM (μ-Capture)

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzer – for the presumptive qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum or plasma (EDTA or heparinized), particularly for women of childbearing age. When performed in conjunction with a Toxoplasma IgG assay, the IMMULITE® 2000 Toxoplasma IgM (μ-Capture) can be used as an aid in the presumptive diagnosis of acute, recent or reactivated *Toxoplasma* infection.

Catalog Number: **L2KTZ2** (200 tests),
L2KTZ6 (600 tests)

Test Code: **TXU** Color: **Violet**

This kit is not FDA cleared for use in the testing (i.e. screening) of blood or plasma donors.

Summary and Explanation

Toxoplasma gondii is an obligate intracellular parasite capable of infecting most mammals, including humans. The organism is transmitted through ingestion of insufficiently cooked meat. Throughout the world, 1 to 90 percent of populations may be infected,⁷ with 25 to 30 percent of the adult population infected in the United States.³ While toxoplasmosis may be manifested in several forms, infections are usually clinically inapparent, and latent infections usually persist for life.² Overt clinical symptoms are similar to infectious mononucleosis, with lymphadenopathy, fever, headache, malaise and sometimes pneumonia and myocarditis.⁶

As with other latent infections, acute toxoplasma infection can pose a serious threat to immunocompromised individuals and newborns who acquire the infection *in utero*. Immunosuppressed patients may develop encephalitis, myocarditis or pneumonitis.³ Congenital infections usually result as a consequence of asymptomatic acute maternal infection.

This infection can cause premature delivery, spontaneous abortion or stillbirth.^{4,6} Neonates may manifest chorioretinitis, hydrocephaly, microcephaly, cerebral calcification and psychomotor retardation.¹ The majority of congenitally infected children will not exhibit any symptoms until later in life.⁴

Management of toxoplasmosis requires serological monitoring of infected individuals,³ as the organism is not readily available for culture. Qualitative testing for the presence of toxoplasma IgM is useful to determine acute infection. Accurate diagnostic information is important, particularly during pregnancy, as treatment with spiramycin can reduce the risk to the fetus.⁵

Principle of the Procedure

IMMULITE 2000 Toxoplasma IgM (μ-Capture) is a solid-phase, two-step chemiluminescent enzyme IgM antibody μ-capture immunoassay. The solid phase, a polystyrene bead, is coated with a monoclonal murine anti-IgM antibody.

The patient specimen and coated bead are added to the Reaction Tube. An alkaline phosphatase-labeled *Toxoplasma* antigen is also added to the Reaction Tube. After the wash and incubation steps, chemiluminescent substrate undergoes hydrolysis in the presence of alkaline phosphatase. IMMULITE 2000 Toxoplasma IgM (μ-Capture) is an immunometric assay. The photon output, as measured by the luminometer, is related to the presence of *Toxoplasma* IgM antibodies in the sample.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The patient need not be fasting, and no special preparations are necessary.

Collect blood aseptically by venipuncture,¹² avoiding hemolysis, into plain, heparinized or EDTA tubes, and separate the serum or plasma from the cells.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Samples which are cloudy or have particulate material should be clarified by low-speed centrifugation.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Toxoplasma IgM (μ -Capture) has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 μ L serum or plasma.

Storage: 3 days at 2–8°C, or 6 months at –20°C.¹⁰

Automatic Predilution Factor: 20.

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Because no test method can offer complete assurance that laboratory specimens do not contain HIV, hepatitis B virus, or other infectious agents, specimens should be handled at the

BSL 2 as recommended for any potentially infectious human serum or blood specimen in the CDC-NIH manual, *Biosafety in Microbiological and Biomedical Laboratories*, 1993.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Toxoplasma IgM (μ -Capture) EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

The Toxoplasma IgM (μ -Capture) Adjustor and Controls supplied with the kit should be fully dissolved after reconstitution. Failure to ensure the homogeneity of the solution may result in poor reproducibility of results.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Toxoplasma IgM (μ -Capture) Bead Pack (L2TZ12)

With barcode. 200 beads, coated with murine monoclonal anti-human IgM antibody. Stable at 2–8°C until expiration date.

L2KTZ2: 1 pack. **L2KTZ6:** 3 packs.

Toxoplasma IgM (μ -Capture) Reagent Wedge (L2TZA2)

With barcode. Two reagents: 11.5 mL of a protein-based buffer. 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to native P-30 *Toxoplasma* antigen, in buffer. Stable at 2–8°C until expiration date.

L2KTZ2: 1 wedge. **L2KTZ6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Toxoplasma IgM (μ -Capture) Adjustor (LTZR)

Lyophilized human serum with IgM reactive to *Toxoplasma*, in buffer, with preservative. The Adjustor serves as the assay's Cutoff. Reconstitute by adding **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. (No further dilution is required.) Stable at 2–8°C for 14 days after reconstitution, or for 6 months (aliquotted) at –20°C.
L2KTZ2: 1 vial. **L2KTZ6:** 2 vials.

Toxoplasma IgM (μ -Capture) Controls (LTZC1, LTZC2)

LTZC1 (Negative Control): One vial containing lyophilized human serum nonreactive to *Toxoplasma*, with preservative. **LTZC2 (Positive Control):** One vial containing lyophilized human serum with IgM reactive to *Toxoplasma*, with preservative. Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 14 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KTZ2: 1 set. **L2KTZ6:** 2 sets.

The IMMULITE 2000 software performs automatic on-board dilution of control samples, and the results will be tracked in the QC database. Enter controls as controls.

For the current control ratio ranges, please refer to the Control insert.

Aliquot Labels with barcodes are supplied with the kit, for use with the Adjustors and Controls. Before use, place the appropriate Aliquot Labels on test tubes, so the barcodes can be read by the barcode reader on the IMMULITE 2000.

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months

(aliquotted) at –20°C.

L2KTZ2: 1 vial. **L2KTZ6:** 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KTZ2: 3 labels **L2KTZ6:** 5 labels

Kit Components Supplied Separately

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2IGZ2: 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2IGZ2: 5 labels.

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also required:

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples:

The Toxoplasma IgM (μ -Capture) Controls (LTZC1-2) supplied with the kit should be used as quality control material to monitor the performance of the assay at the cutoff range. The Positive Control is used to

validate the IMMULITE 2000 Toxoplasma IgM (μ-Capture) assay at a critical level when determining the presence of an active toxoplasma infection.

In addition to the controls provided, users may wish to run additional controls for their own purposes

Additional controls may be tested in accordance with guidelines or requirements of local, state, and/or federal regulations or accrediting organizations.

It is also recommended that known reactive and nonreactive specimens be run periodically to assure pipetting accuracy for the dilution step.

The controls should be processed at (or near) the beginning of every run containing patient samples to be tested for *Toxoplasma* IgM, and also when readjusting.

Users of the IMMULITE 2000 system are advised to consult the NCCLS document C24-A, Internal Quality Control Testing: Principles and Definitions, for additional guidance on quality control for the basic principles and definitions in dealing with internal quality control testing.

Calculation of Cutoff and S/CO Ratio:
The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times \text{P1}}$$

Calculation and reporting of qualitative (reactive/nonreactive/indeterminate) and s/co ratio results are handled automatically by the IMMULITE 2000.

The result for a sample is reported as "Indeterminate" if the counts per second for that sample fall within $\pm 10\%$ of the cutoff. The result is reported as "Reactive"

if the sample's counts are above the indeterminate range, and "Nonreactive" if below this range.

Interpretation of Results

The cutoff of the IMMULITE 2000 Toxoplasma IgM (μ-Capture) assay was determined with reactive and nonreactive patient samples by a ROC analysis with a balanced consideration of sensitivity and specificity.

A result of **"Reactive"** (ratio of ≥ 1.1) indicates that the patient's sample is reactive and that *Toxoplasma* IgM antibodies were detected in the sample.

A result of **"Nonreactive"** (ratio of < 0.9) indicates that the patient's sample is nonreactive and that *Toxoplasma* IgM antibodies were not detected in the sample.

Any result of **"Indeterminate"** (ratio between 0.9 and < 1.1) should be retested. Samples which still test as "Indeterminate" should be tested by an alternate method, or a second sample should be taken — if possible — within a reasonable period of time (e.g., one week).

The presence of IgM antibodies to *Toxoplasma gondii* is an indication of a recent exposure to the organism.¹¹

The magnitude of the measured results (cps) above the Cutoff is not indicative of the total amount of antibodies detected.

An interpretation of the results is not recommended if the positive or negative control falls outside the range specified in the Quality Control section.

Reports by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Toxoplasma IgM (μ-Capture) EIA. Values obtained from other manufacturers' assay methods may not be used interchangeably."

Diagnosis of acute infection should not rely on any single test result. If acute infection is suspected, a patient sample should be tested for the presence of both *Toxoplasma*-specific IgG and IgM antibodies. The results should be interpreted according to the guidelines provided in the following table and a new sample obtained, if applicable.

Anti- <i>Toxoplasma gondii</i> Result		Report/Interpretation
IgM	IgG	
Nonreactive	Nonreactive	It is presumed the patient has not been infected with and is not undergoing an acute infection with <i>Toxoplasma gondii</i> . If symptoms persist submit a new specimen within three weeks.
Nonreactive	Reactive	From this testing it cannot be determined whether the patient is or is not undergoing a reactivated <i>Toxoplasma gondii</i> infection. It appears the patient has been previously infected with <i>Toxoplasma gondii</i> . Infection occurred more than one year ago.
Nonreactive	Equivocal	Obtain a new specimen for further testing. The patient may not be undergoing an acute infection with <i>Toxoplasma gondii</i> . Determining whether the patient has been previously infected with <i>Toxoplasma gondii</i> is not possible.
Equivocal	Nonreactive	Obtain a new specimen for determination of IgM antibodies to <i>Toxoplasma gondii</i> . It cannot be determined if the patient is undergoing an acute <i>Toxoplasma gondii</i> infection. It appears the patient has not been previously infected with <i>Toxoplasma gondii</i> . If the new specimen result is reactive or equivocal for IgM antibodies, the specimen should be sent to a reference laboratory with experience in the diagnosis of toxoplasmosis for further testing.
Equivocal	Reactive	Obtain a new specimen for determination of IgM antibodies to <i>Toxoplasma gondii</i> . It cannot be determined if the patient is undergoing or has undergone an acute <i>Toxoplasma gondii</i> infection. It appears the patient has been previously infected with <i>Toxoplasma gondii</i> . If the new specimen result is equivocal or reactive for IgM antibodies, the specimen should be sent to a reference laboratory with experience in the diagnosis of toxoplasmosis for further testing.
Equivocal	Equivocal	Obtain a new specimen for further testing. It cannot be determined if the patient is undergoing an acute infection or has been previously infected with <i>Toxoplasma gondii</i> . If the new specimen result is equivocal or reactive for IgM antibodies, the specimen should be sent to a reference laboratory with experience in the diagnosis of toxoplasmosis for further testing.

Anti- <i>Toxoplasma gondii</i> Result		Report/Interpretation
IgM	IgG	
Reactive	Nonreactive	Obtain a new specimen for further testing. The patient may or may not be acutely infected with <i>Toxoplasma gondii</i> . Since the IgG antibodies to <i>Toxoplasma gondii</i> are nonreactive, the specimen may have been obtained too early in the disease process for an accurate determination. Retest the new specimen with a different anti- <i>Toxoplasma gondii</i> IgM assay. If the new specimen result is still reactive for IgM antibodies, the specimen should be sent to a reference laboratory with experience in the diagnosis of toxoplasmosis for further testing.
Reactive	Reactive	The patient may or may not be acutely infected with <i>Toxoplasma gondii</i> . Obtain a new specimen for further testing. Since the IgG antibodies to <i>Toxoplasma gondii</i> are reactive, it appears the patient may be acutely infected with <i>Toxoplasma gondii</i> . The new specimen should be repeated with a different anti- <i>Toxoplasma gondii</i> IgM assay. If the new specimen result is still reactive for IgM and IgG antibodies to <i>Toxoplasma gondii</i> , the specimen should be sent to a reference laboratory with experience in the diagnosis of toxoplasmosis for further testing.
Reactive	Equivocal	It cannot be determined if the patient is acutely infected with <i>Toxoplasma gondii</i> . Obtain a new specimen for further testing. Determining whether the patient has been previously infected with <i>Toxoplasma gondii</i> is not possible. The specimen may have been collected too early during the disease process for an accurate determination. Retest the new specimen with a different anti- <i>Toxoplasma gondii</i> IgM assay. If the new specimen result is still reactive for IgM and the IgG is reactive/nonreactive/equivocal for antibodies to <i>Toxoplasma gondii</i> , the specimen should be sent to a reference laboratory with experience in the diagnosis of toxoplasmosis for further testing.

Expected Values

Individuals infected with the *Toxoplasma* organism will typically exhibit detectable levels of IgM antibody immediately before or soon after the onset of symptoms.² IgM titers normally decline within four to six months, but may persist at low levels up to a year.⁴ Patients with active *toxoplasma* chorioretinitis usually have undetectable levels of IgM.⁴

The prevalence of *Toxoplasma* infection can vary depending on a number of factors such as age, gender, geographical location, socio-economic status, race, type of test used, specimen collection and handling procedures, and clinical and epidemiological history of individual patients. There are approximately 3,000 cases of congenital toxoplasmosis reported per year, with an average of 0.6 cases per 1,000 pregnancies in the United States.¹³

Limitations

The results of the test must be taken within the context of the patient's clinical history, symptomology and other laboratory findings. Results of the IMMULITE 2000 *Toxoplasma* IgM (μ-Capture) kit are not by themselves diagnostic and should be interpreted in light of the patients' clinical condition and the results of other diagnostic procedures.

For the determination of seroconversion from *nonreactive* to *reactive*, two serum samples should be drawn three to four weeks apart, during the acute and convalescent stages of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution. A nonreactive result for *Toxoplasma* IgM does not preclude the possibility of an acute infection in immunocompromised patients. *Toxoplasma gondii*-specific IgG antibodies are generally low, and *Toxoplasma gondii*-specific IgM antibodies may be undetectable in patients who are immunocompromised.

The performance characteristics of this assay have not been established for use with neonates, cord blood, or pretransplant patients.

Antibodies to human IgG have been added to the reagents to remove specific IgG and rheumatoid factors which may cause false reactive responses.

The performance of the IMMULITE 2000 *Toxoplasma* IgM (μ-Capture) assay has not been established for sample matrices other than human serum and plasma.

The presence of IgM antibodies in a single specimen is not sufficient to distinguish between active or past infection. Patients suspected of having primary or active infection should be tested for the presence of IgG antibodies to *Toxoplasma gondii*.

If a treatment is prescribed early enough, antibody production is decreased, and IgG and IgM levels remain low and can coexist for years.

The continued presence or absence of antibodies cannot be used to determine the success or failure of therapy.

Reactive test results may not be valid in persons who have received blood transfusions or other blood products within the past several months.

Testing should not be performed as a screening procedure for the general population. The predictive value of a reactive or nonreactive serologic result depends on the pre-test likelihood of toxoplasmosis being present. Testing should only be done when clinical evidence suggests the diagnosis of toxoplasmosis.

This test is not intended for the determination of immune status. It is intended for the determination of a patient's antibody response to indicate active infection to *Toxoplasma gondii* and not as an indication of immunity.

Low levels of IgM antibodies may occasionally persist for more than 12 months post-infection. Such a residual antibody response may be distinguished from the early IgM response to active infection by testing sera from the patient 2-4 weeks later using the IMMULITE 2000 *Toxoplasma* IgM (μ-Capture) kit and by reference to changing levels of *Toxoplasma* IgG antibodies.

Specific IgM antibodies are usually detected in patients with recent primary infection, but they may be found in patients with reactivated or secondary infections, and they are sometimes found in patients with no other detectable evidence of recent infection.

Samples taken too early during the course of a primary infection with *Toxoplasma gondii* may not contain detectable levels of IgM-specific antibody. In some patients, IgM-specific antibody results may revert to nonreactive levels within three weeks after infection with *Toxoplasma gondii*.

Measurement of *Toxoplasma gondii*-specific IgG antibodies may also be of some value in the serological assessment of these patients.

With very low prevalence analytes, such as anti-*Toxoplasma gondii* IgM, there is the increased possibility that a reactive result is actually a false reactive, reducing the assay's positive predictive value.

Due to the high sensitivity of the assay, samples from a *Toxoplasma* subacute infection may be detectable. Results from these samples should be further evaluated in the context of the clinical examination, patient medical history, and other findings.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed as a signal-to-cutoff ratio.

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Crossreactivity: A study was conducted to evaluate whether the measurement of *Toxoplasma* IgM antibody is affected by closely related microorganisms. Eighty-seven seronegative sera containing antibodies to varicella zoster virus ($n=3$), measles ($n=10$), cytomegalovirus (CMV) ($n=10$), herpes simplex virus ($n=10$), *Toxoplasma* ($n=10$), *Mycoplasma pneumoniae* ($n=10$), Epstein-Barr virus ($n=10$), syphilis ($n=10$), parvovirus ($n=8$) and rheumatoid factor ($n=6$) were tested by IMMULITE 2000 *Toxoplasma* IgM (μ -Capture), and all yielded nonreactive results.

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 539 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL may have an effect on results. (See "Lipemia" table.)

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 34 volunteers into plain, heparinized, EDTA and SST[®] Becton Dickinson plastic vacutainer tubes. Eleven samples were spiked with sera containing IgM antibodies to *Toxoplasma*, and all samples were assayed by the IMMULITE 2000 *Toxoplasma* IgM (μ -Capture) procedure. Results are expressed as a signal-to-cutoff ratio. By linear regression:

(EDTA Plastic) = 1.02 (Serum Plastic) + 0.009
 $r = 0.998$

(Heparin Plastic) = 1.00 (Serum Plastic) – 0.002
 $r = 0.996$

(SST Plastic) = 1.03 (Serum Plastic) – 0.006
 $r = 0.994$

Means:

0.53 (Serum Plastic)

0.55 (EDTA Plastic)

0.54 (Heparin Plastic)

0.54 (SST Plastic)

In another experiment, blood was collected from 40 volunteers into glass and plastic serum Becton Dickinson vacutainer tubes. All samples were assayed by the IMMULITE 2000 Toxoplasma IgM (μ -Capture) procedure. Results are expressed as a signal-to-cutoff ratio. By linear regression:

$$(\text{Serum Plastic}) = 1.00 (\text{Serum Glass}) + 0.010$$
$$r = 0.808$$

Means:

0.11 (Serum Plastic)

0.10 (Serum Glass)

Method Comparison: The assay was compared to a commercially available μ -capture assay for *Toxoplasma* IgM (Kit A) on 195 samples:

Kit A	IMMULITE 2000 Toxoplasma IgM (μ -Capture) (L2KKTZ)		
	Pos	Ind	Neg
Reactive	18	1	0
Indeterminate	1	0	1
Nonreactive	4	0	170

Total Agreement: 96.4% (188/195)

Relative Sensitivity: 94.7% (18/19)

Relative Specificity: 97.7% (170/174)

Positive Predictive Value: 78.3% (18/23)

Negative Predictive Value: 99.4% (170/171)

Clinical Performance

In a clinical study performed in Europe, specimens from 452 defined *Toxoplasma* IgM negative patients and 40 defined *Toxoplasma* acute phase infection patients were evaluated by the IMMULITE 2000 Toxoplasma IgM (μ -Capture) procedure. The specimens were analyzed and *Toxoplasma* IgM positive and negative indications were defined by well established diagnostic methods, clinical findings and serological patterns. Of the 452 patients defined as *Toxoplasma* IgM negative, 396 were positive for *Toxoplasma* IgG, reflecting a latent or distant infection. Among these 396 latent infection patients, 293 were collected from pregnant women representing all stages of pregnancy. The tables below present the results of this study.

Comparison for All Subjects:

Clinically Defined Samples

IMMULITE 2000 Toxoplasma IgM μ -Capture (L2KKTZ)	Clinically Defined Samples		
	Pos	Ind	Neg
Reactive	40	0	20
Indeterminate	0	0	6
Nonreactive	0	0	426

With Indeterminate Cases:

Agreement: 95.9% (472/492); Sensitivity: 100% (40/40); Specificity: 95.6% (432/452)

Without Indeterminate Cases:

Agreement: 95.9% (466/486); Sensitivity: 100% (40/40); Specificity: 95.5% (426/446)

Comparison for Pregnant Subjects:

Clinically Defined Samples

IMMULITE 2000 Toxoplasma IgM μ -Capture (L2KKTZ)	Clinically Defined Samples		
	Pos	Ind	Neg
Reactive	0	0	14
Indeterminate	0	0	5
Nonreactive	0	0	274

With Indeterminate Cases:

Specificity: 95.2% (279/293)

Without Indeterminate Cases:

Specificity: 95.1% (274/288)

References

- 1) Daffos F. Prenatal management of 746 pregnancies at risk for congenital toxoplasmosis. N Engl J Med 1988;318: 271-5.
- 2) Krick JA, Remington JS. Toxoplasmosis in the adult: an overview. N Engl J Med 1978;298: 550-3.
- 3) Krogstad DJ, et al. Blood and tissue protozoa. In: Lennette EH, et al, editors. Manual of clinical microbiology. 4th ed. Washington, D.C.: American Society for Microbiology, 1985: 612-30.
- 4) McCabe RE, Remington JS. *Toxoplasma gondii*. In: Mandell, Douglas, Bennett, editors. Principles and practice of infectious diseases. 2nd ed. New York: John Wiley, 1985.
- 5) Jeannel D, et al. What is known about the prevention of congenital toxoplasmosis? Lancet 1990;336:359-61.
- 6) Remington JS, McLeod R. Toxoplasmosis. In: Braude, editor. Microbiology and infectious diseases. Philadelphia: Saunders, 1981.
- 7) Walls KW. Toxoplasmosis. In: Balows A, et al, editors. Laboratory diagnosis of infectious diseases. New York: Springer-verlag, 1988: 998-1017.
- 8) Baltz ML, Searcy RL. Clinical significance and advanced serologic diagnosis of ToRCH infections. Am Clin Lab 1994;March/April:18-23.
- 9) Wilson M, et al. Seroepidemiology of toxoplasmosis in the

United States. [Reprint from Abstracts presented at the 1994 American Society for Microbiology, May 1994]. 10) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:354-60. 11) van Loon AM, van der Logt JTM, et al. Enzyme-linked immunosorbent assay that uses labeled antigen for detection of immunoglobulin M and A antibodies in Toxoplasmosis: comparison with indirect immunofluorescence and double-sandwich enzyme-linked immunosorbent assay. J Clin Microbiol 1983; June: 997-1004. 12) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 13) Isada, CM, et al. Infectious Diseases Handbook. Lexi-Comp and American Pharmaceutical Association. Hudson, OH. 4th Edition, 2001:317

Technical Assistance

For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.
©2009 Siemens Healthcare Diagnostics Inc. All rights reserved.