NT-proBNP

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 NT-proBNP

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the quantitative measurement of N-terminal pro-brain natriuretic peptide (NT-proBNP) in heparinized plasma, as an aid in the diagnosis of individuals suspected of having congestive heart failure. The test is further indicated for the risk stratification of patients with acute coronary syndrome and congestive heart failure.

Catalog Numbers: **L2KNT2** (200 tests)

Test Code: **NTP** Color: **Red**

Summary and Explanation

Chronic heart failure, most commonly caused by left ventricular dysfunction, is a clinical syndrome resulting from reduced heart pumping capacity. It currently affects about 5 million Americans and is the leading cause of hospitalization for patients older than age 65.¹

Studies have shown that B-type natriuretic peptide (BNP), a member of the natriuretic peptide family, is useful in assisting the diagnosis of heart failure.² ProBNP, the 108 amino acid precursor of BNP, is synthesized and released by ventricular myocytes in response to increased ventricular wall tension. ProBNP is cleaved to release a 32 amino acid active hormone (BNP), and the N-terminal 76 amino acid inactive peptide NT-proBNP. BNP antagonizes the renin-angiotensin-aldosterone (RAAS) system by causing natriuresis, diuresis and vasodilation — inhibiting renin excretion and aldosterone production.³ The plasma levels of both BNP and NT-proBNP are markedly increased in subjects with left ventricular dysfunction, and correlate well with the New York Heart Association functional classification of heart failure.^{5,6}

NT-proBNP has shown advantages over BNP as a biochemical marker because of its longer half-life, better *in vitro* stability, reduced intra-individual fluctuation, and higher circulating concentration.⁴ Determination of NT-proBNP levels in

plasma can be useful in identifying patients with chronic heart failure, assessing the severity, predicting increased morbidity, and monitoring the therapeutic response.^{1,4,7}

Principle of the Procedure

IMMULITE 2000 NT-proBNP is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

Serum is not suitable for use.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 NT-proBNP has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL heparinized plasma.

Storage: 3 days at 2–8°C or 6 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
	Contains: sodium azide; NT-proBNP Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

NT-proBNP Bead Pack (L2NT12)

With barcode. 200 beads, coated with polyclonal sheep anti-NT-proBNP. Stable at 2–8°C until expiration date.

L2KNT2: 1 pack.

NT-proBNP Reagent Wedge (L2NTA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal sheep anti-NT-proBNP in buffer. Stable at 2–8°C until expiration date.

L2KNT2: 1 wedge.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

NT-proBNP Adjustors (LNTL, LNTH)

Two vials (Low and High) of lyophilized NT-proBNP in a nonhuman serum matrix, with preservative. Reconstitute each vial by adding **3.0 mL** distilled or deionized water. Mix by gentle swirling or inversion. Stable at –20°C for 3 months (aliquotted) after reconstitution.

L2KNT2: 1 set.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

NT-proBNP Sample Diluent (L2NTZ)

For the on-board dilution of patient samples. One vial of concentrated (ready-to-use) processed, NT-proBNP-free nonhuman serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2NTZ: 25 mL.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2NTZ: 3 labels.

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LNTCM: A bi-level, nonhuman serum-based NT-proBNP Control Module.

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of NT-proBNP.

Expected Values

The reference range values were generated on the IMMULITE 2500 analyzer. Based on its relationship to IMMULITE 2500 NT-proBNP (See Method Comparison 2), IMMULITE 2000 NT-proBNP can be expected to have essentially the same reference range.

A reference range study was performed in which 217 heparinized plasma samples — 124 from apparently healthy laboratory volunteers and 93 from presumed normal specimens from a large reference laboratory, both with no known cardiac disease history — were assayed by the IMMULITE 2500 NT-proBNP procedure. Medians and 95th percentiles for relevant age groups are tabulated below. The most appropriate decision thresholds are 125 pg/mL for patients younger than 75 years old, and 450 pg/mL for patients 75 years or older.⁴

Reference Group		
	<75 yrs (n = 124)	75+ yrs (n = 93)
Median	28.5 pg/mL	172 pg/mL
95 th Percentile	110 pg/mL	589 pg/mL

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Heterophilic antibodies in human serum/plasma can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on heparinized plasma samples collected in tubes without gel barriers.)

Conversion Factor:

pg/mL \times 0.118 \rightarrow pmol/L

Reportable Range: 20 – 35,000 pg/mL
(2.4 – 4,130 pmol/L).

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 10 pg/mL
(1.2 pmol/L).

High-dose Hook Effect:

None up to 385,000 pg/mL.

Precision: Samples were assayed in duplicate over the course of 10 days, four runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: It is recommended that only samples above the assay's upper limit (> 35,000 pg/mL or > 4,130 pmol/L) should be diluted. Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three NT-proBNP solutions (26,500, 53,000, and 106,000 pg/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for NT-proBNP. (See "Specificity" table.)

Bilirubin: May interfere with the assay in NT-proBNP samples below 500 pg/mL, causing degradation of values. (See "Bilirubin" table for representative data.)

Hemolysis: Presence of hemoglobin in concentrations up to 600 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 25 volunteers into plastic lithium-heparin and plastic plasma gel barrier tubes (PST®). All tubes were

from Becton Dickinson. In 17 sets of the matched samples, equal volumes were spiked with various concentrations of NT-proBNP, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 NT-proBNP procedure.

(PST[®] plasma) = 0.991 (Li-Heparin plasma) + 138 pg/mL
r = 0.994

Means:

5,457 pg/mL (Li-Heparin plasma)
5,548 pg/mL (PST[®] plasma)

Method Comparison 1: The assay was compared to a commercially available Elecsys 2010 proBNP kit on 366 patient samples. (Concentration range: approximately 7.9 to 34,468 pg/mL. See graph 1.) By linear regression:

(IML 2000) = 0.97 (Elecsys 2010) – 9.9 pg/mL
r = 0.992

Means:

1,978 pg/mL (IMMULITE 2000)
2,029 pg/mL (Elecsys 2010)

Method Comparison 2: The assay was compared to IMMULITE 2500 NT-proBNP on 366 patient samples. (Concentration range: 21.3 to 32,855 pg/mL. See graph 2.) By linear regression:

(IML 2000) = 1.01 (IML 2500) – 35.8 pg/mL
r = 0.996

Means:

2,231 pg/mL (IMMULITE 2000)
2,246 pg/mL (IMMULITE 2500)

References

1) Cardarelli R, et al. B-type Natriuretic Peptide: a review of its diagnostic, prognostic, and therapeutic monitoring value in heart failure for primary care physicians. J Am Board Fam Pract. 2003;16(4):327-33. 2) Richarde AM, et al. Plasma N-terminal pro-Brain natriuretic peptide and adrenomedullin: new neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998;97:1921-29. 3) De Lemos JA, et al. B-type natriuretic peptide in cardiovascular disease. The Lancet 2003;362:316-22. 4) McCullough PA, et al. Sorting out the evidence on natriuretic peptides. Rev Cardiovasc Med. 2003;4:S13-S19. 5) Hunt PJ, et al. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): a new marker for cardiac impairment. Clin Endocrinol 1997;47:287-96. 6) Talwar S, et al. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left-ventricular systolic dysfunction in a high risk population. Eur Heart J 1999;20:1736-44. 7) Gardner RS, et al. N-terminal pro-brain

natriuretic peptide: a new gold standard in predicting mortality in patients with advanced heart failure. Eur Heart J 2003;24:1735-43.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

Manufactured under license from Roche Diagnostics GmbH.

Tables and Graphs

Precision (pg/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	35.6	1.92	5.4%	2.29	6.4%
2	146	4.58	3.1%	6.46	4.4%
3	519	16.6	3.2%	23.9	4.6%
4	1,430	42.4	3.0%	57.0	4.0%
5	8,884	204	2.3%	435	4.9%
6	15,711	375	2.4%	643	4.1%
7	29,725	1,211	4.1%	1,400	4.7%

Dansk. ¹I samme kørsel, ²Total, ³Middelværdi, ⁴Standardafvigelse, ⁵Varianskoefficient. **Eesti.** ¹Mõõtmisseeria sisene, ²Seeriaveheline, ³Keskmine väärtused, ⁴Standardhälve, ⁵Variatsioonikoeffitsient. **Latviski.** ¹Sērijas robežās, ²Kopā, ³Vidējā vērtība, ⁴SD, ⁵Variāciju koeficients. **Lietuviškai.** ¹Vieno tyrimo metu, ²Bendrai, ³Vidurkiai, ⁴SD, ⁵CV. **Norsk.** ¹Innen serien, ²Total, ³Middel, ⁴SD, ⁵CV. **Svenska.** ¹Inom körning, ²Totalt, ³Medel, ⁴SD, ⁵CV.

Linearity (pg/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	31,015	—	—
	4 in 8	18,327	15,508	118%
	2 in 8	9,418	7,754	121%
	1 in 8	4,326	3,877	112%
2	8 in 8	30,139	—	—
	4 in 8	17,650	15,070	117%
	2 in 8	8,769	7,535	116%
	1 in 8	4,103	3,767	109%
3	8 in 8	28,681	—	—
	4 in 8	16,041	14,340	112%
	2 in 8	7,848	7,170	109%
	1 in 8	3,566	3,585	99%
4	8 in 8	26,741	—	—
	4 in 8	15,628	13,370	117%
	2 in 8	7,861	6,685	118%
	1 in 8	3,547	3,343	106%
5	8 in 8	25,288	—	—
	4 in 8	13,787	12,644	109%
	2 in 8	6,503	6,322	103%
	1 in 8	3,063	3,161	97%
6	8 in 8	24,542	—	—
	4 in 8	14,246	12,271	116%
	2 in 8	6,956	6,135	113%
	1 in 8	3,169	3,068	103%

Dansk. ¹Fortynding, ²Observeret, ³Forventet, ⁴%Obs./ Forv., ⁵Ufortyndet. **Eesti.** ¹Lahjendus, ²Täheldatud (T), ³Oodatud (O), ⁴% T/O, ⁵8 8-st. **Latviski.** ¹Atšķaidījums, ²Novērots (N), ³Gaidīts (G), ⁴% N/G, ⁵8 : 8. **Lietuviškai.** ¹Skiedimas, ²Nustatyta (N), ³Tikėtasi (T), ⁴% N/T, ⁵8 iš 8. **Norsk.** ¹Fortynning, ²Observert (O), ³Forventet (E), ⁴% O/E, ⁵8 in 8. **Svenska.** ¹Spädning, ²Observerat (O), ³Förväntat (E), ⁴% O/E, ⁵Ospätt.

Recovery (pg/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	42.0	—	—
	A	1,256	1,365	92%
	B	2,474	2,690	92%
	C	5,036	5,340	94%
2	—	634	—	—
	A	1,972	1,928	102%
	B	3,414	3,253	105%
	C	6,016	5,903	102%
3	—	2,147	—	—
	A	3,460	3,365	103%
	B	4,932	4,690	105%
	C	7,839	7,340	107%
4	—	4,102	—	—
	A	6,845	5,222	131%
	B	8,402	6,547	128%
	C	8,848	9,197	96%
5	—	11,749	—	—
	A	12,406	12,486	99%
	B	14,185	13,811	103%
	C	16,302	16,461	99%
6	—	15,521	—	—
	A	15,913	16,070	99%
	B	16,403	17,395	94%
	C	19,744	20,045	98%

Dansk. ¹Opløsning, ²Observeret, ³Forventet, ⁴%Obs./ Forv. **Eesti.** ¹Lahus, ²Täheldatud (T), ³Oodatud (O), ⁴% T/O. **Latviski.** ¹Šķīdums, ²Novērots (N), ³Gaidīts (G), ⁴% N/G. **Lietuviškai.** ¹Tirpalas, ²Nustatyta (N), ³Tikėtasi (T), ⁴% N/T. **Norsk.** ¹Løsning, ²Observert (O), ³Forventet (E), ⁴% O/E. **Svenska.** ¹Lösning, ²Observerat (O), ³Förväntat (E), ⁴% O/E.

Specificity

Compound ¹	Amount Added ² (pg/mL)	Apparent ³ pg/mL	% Cross reactivity ⁴
ANP	3,100,000	<20 pg/mL	ND
NT-proANP (1-30)	3,500,000	<20 pg/mL	ND
NT-proANP (31-67)	1,000	<20 pg/mL	ND
NT-proANP (79-98)	1,000	<20 pg/mL	ND
BNP32	3,500,000	<20 pg/mL	ND
CNP32	2,200,000	<20 pg/mL	ND
Adrenomedullin	1,000	<20 pg/mL	ND
Aldosterone	600	<20 pg/mL	ND
Angiotensin I	600	<20 pg/mL	ND
Angiotensin II	600	<20 pg/mL	ND
Angiotensin III	1,000	<20 pg/mL	ND
Endothelin	20	<20 pg/mL	ND
Renin	50,000	<20 pg/mL	ND
Urodilatin	3,500,000	<20 pg/mL	ND
Arg-Vasopressin	1,000	<20 pg/mL	ND

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

IMMULITE is a trademark of
Siemens Healthcare Diagnostics.

©2008 Siemens Healthcare Diagnostics.
All rights reserved.

Made in: UK

Siemens Healthcare
Diagnostics Products Ltd.
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

2015-06-17

PINL2KNT – 2 {16}

Changes in this Edition:

cc#EU22565: 1) Removed references to discontinued products. 2) Under Warnings and Precautions, added caution statement, required GHS pictograms/signal words/hazard and precautionary codes and statements, offending chemicals and sources. 3) Changed "Origin: UK" to "Made in: UK". 4) Updated the manufacturer address to include Glyn Rhonwy. 5) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

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