PAPP-A

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 PAPP-A

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the quantitative measurement of pregnancy-associated plasma protein A (PAPP-A) in serum or heparinized plasma.

Catalog Number: **L2KPC2** (200 tests)

Test Code: **PAA** Color: **Brown**

Summary and Explanation

Pregnancy-associated plasma protein (PAPP-A) is a large molecular weight placenta-derived glycoprotein.¹ During pregnancy it is produced in high concentrations by the trophoblast and released into maternal circulation. PAPP-A levels in maternal serum rise steadily with gestational age, most noticeably during the last part of pregnancy.² The functional significance of PAPP-A is unclear. However, some studies suggest that reduced PAPP-A concentrations in pregnancy are associated with chromosomal abnormalities in the fetus.³

Maternal serum PAPP-A assessment, between 11 and 14 weeks of pregnancy, is reported to have significant utility in screening for Down syndrome and other chromosomal anomalies. A combination of maternal age-related risk, free β -HCG and fetal nuchal translucency (NT) measurements might substantially increase the efficiency of prenatal screening compared to second-trimester screening; using this approach, various investigators have reported detection rates for Down syndrome of 85 to 90 percent at a 5 percent false-positive rate.^{6,7}

Principle of the Procedure

IMMULITE 2000 PAPP-A is a solid-phase, enzyme-labeled chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

EDTA plasma is not suitable for use.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 PAPP-A has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 μ L serum or plasma.

Storage: 24 hours at 2–8°C or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁸⁻¹⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: 2-methyl-2H-isothiazol-3-one; PAPP-A Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

PAPP-A Bead Pack (L2PC12)

With barcode. 200 beads, coated with monoclonal murine anti-PAPP-A. Stable at 2–8°C until expiration date.

L2KPC2: 1 pack

PAPP-A Reagent Wedge (L2PCA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-PAPP-A in buffer. Stable at 2–8°C until expiration date.

L2KPC2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

PAPP-A Adjustors (LPCL, LPCH)

Two vials (Low and High) of lyophilized PAPP-A, in a nonhuman serum matrix. Reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 30 days after reconstitution, or for 3 months (aliquotted) at –20°C.

L2KPC2: 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 2 (L2M2Z, L2M2Z4)

For the on-board dilution of patient samples. One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes
(16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LPCCM: Bi-level control module in a nonhuman serum matrix.

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 4 weeks.

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of PAPP-A.

Expected Values

Each laboratory should establish its own reference ranges.

Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: Up to 10 mIU/mL

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 0.025 mIU/mL

High-Dose Hook Effect: None up to 115 mIU/mL.

Precision: Samples were repeatedly assayed in duplicate over the course of several days, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three PAPP-A solutions (20, 40 and 80 mIU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for PAPP-A. (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Effects of the presence of hemoglobin were evaluated by spiking 10, 20, and 30 μ L of a red cell hemolysate (corresponding to 157, 314, and 550 mg/dL hemoglobin) into 1 mL aliquots of 5 patient samples. Hemoglobin in concentrations up to 157 mg/dL has no effect on results, within the precision of the assay. Higher levels of hemoglobin may cause slight elevations in results (See "Hemoglobin" table).

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain and heparinized tubes. Equal volumes of the matched samples were spiked with various concentrations of PAPP-A, to obtain values throughout the calibration

range of the assay, and then assayed by the IMMULITE 2000 PAPP-A procedure.

(Heparin) = 1.08 (Serum) + 0.08 mIU/mL
r = 0.996

Means:

4.72 mIU/mL (Serum)

5.18 mIU/mL (Heparin)

In another study, blood was collected from 24 volunteers into plain and EDTA tubes. Equal volumes of the matched samples were spiked with various concentrations of PAPP-A, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 PAPP-A procedure.

(EDTA) = 1.21 (Serum) + 0.02 mIU/mL
r = 0.995

Means:

3.65 mIU/mL (Serum)

4.43 mIU/mL (EDTA)

EDTA plasma is not suitable for use.

Method Comparison 1: The assay was compared to a commercially available assay for PAPP-A (Kit A) on 152 patient samples. (Concentration range: approximately 0.3 to 10 mIU/mL. See graph 1.) By linear regression:

(IML 2000) = 0.99 (Kit A) + 0.5 mIU/mL
r = 0.973

Means:

3.82 mIU/mL (IMMULITE 2000)

3.36 mIU/mL (Kit A)

Method Comparison 2: The assay was compared to IMMULITE PAPP-A on 283 samples. (Concentration range: approximately 0.3 to 10 mIU/mL. See graph 2.) By linear regression:

(IML 2000) = 1.03 (IML) – 0.06 mIU/mL
r = 0.995

Means:

3.70 mIU/mL (IMMULITE 2000)

3.64 mIU/mL (IMMULITE)

References

1) Lin, T-M, *et al.* Characterization of four human pregnancy-associated plasma proteins. *Am J Obstet Gynecol* 1974; 118: 223-236. 2) Lin, T-M, *et al.* Measurement of pregnancy-associated plasma proteins during human gestation. *J Clin Invest* 1974; 54: 576-82. 3) Pederson, J. F. *et al.* Serum levels of human placental lactogen, pregnancy-associated plasma protein A and endometrial secretory PP14 in first trimester of diabetic pregnancy. *Acta Obstet Gynecol Scand* 1998; 77:155-8. 4) Wald NJ, *et al.* Serum screening for Down's syndrome between 8 and

14 weeks of pregnancy. *Br J Obstet Gynaecol* 1996; 103: 407-12. 5) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 6) Bindra R, *et al.* One-stop clinic for assessment of risk for trisomy 21 at 11-14 weeks: a prospective study of 15,030 pregnancies. *Ultrasound Obstet Gynecol* 2002;3:219-25. 7) Wald NJ, Hackshaw AK. Advances in antenatal screening for Down syndrome. *Baillieres Best Pract Res Clin Obstet Gynaecol* 2000;14(4):563-80. 8) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8. 9) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 10) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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Linearity (mIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	1.25	—	—
	4 in 8	0.62	0.63	98%
	2 in 8	0.29	0.31	94%
	1 in 8	0.14	0.16	88%
2	8 in 8	2.49	—	—
	4 in 8	1.22	1.25	98%
	2 in 8	0.63	0.62	102%
	1 in 8	0.29	0.31	94%
3	8 in 8	4.43	—	—
	4 in 8	2.25	2.22	101%
	2 in 8	1.08	1.11	97%
	1 in 8	0.55	0.55	100%
4	8 in 8	5.08	—	—
	4 in 8	2.62	2.54	103%
	2 in 8	1.32	1.27	104%
	1 in 8	0.62	0.64	97%
5	8 in 8	8.17	—	—
	4 in 8	4.27	4.09	104%
	2 in 8	2.23	2.04	109%
	1 in 8	1.08	1.02	106%

Dansk. ¹Fortynding, ²Observeret, ³Forventet, ⁴% Obs./ Forv., ⁵Ufortyndet. **Eesti.** ¹Lahjendus, ²Täheldatud (T), ³Oodatud (O), ⁴% T/O, ⁵8 8-st. **Latviski.** ¹Atšķaidījums, ²Novērots (N), ³Gaidīts (G), ⁴% N/G, ⁵8 : 8. **Lietuviškai.** ¹Skiedimas, ²Nustatyta (N), ³Tikėtasi (T), ⁴% N/T, ⁵8 iš 8. **Svenska.** ¹Spädning, ²Observerat (O), ³Förväntat (E), ⁴% O/E, ⁵Ospätt.

Recovery (mIU/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	ND	—	—
	A	1.38	1.0	138%
	B	2.26	2.0	113%
	C	4.01	4.0	100%
2	—	1.1	—	—
	A	2.32	2.05	113%
	B	3.28	3.05	108%
	C	5.52	5.05	109%
3	—	3.43	—	—
	A	4.30	4.26	101%
	B	5.48	5.26	104%
	C	7.20	7.26	99%
4	—	4.47	—	—
	A	5.24	5.25	100%
	B	6.21	6.25	99%
	C	8.21	8.25	100%
5	—	6.70	—	—
	A	7.21	7.37	98%
	B	8.21	8.37	98%
	C	9.95	10.37	96%

Dansk. ¹Opløsning, ²Observeret, ³Forventet, ⁴% Obs./ Forv. **Eesti.** ¹Lahus, ²Täheldatud (T), ³Oodatud (O), ⁴% T/O. **Latviski.** ¹Šķīdums, ²Novērots (N), ³Gaidīts (G), ⁴% N/G. **Lietuviškai.** ¹Tirpalas, ²Nustatyta (N), ³Tikėtasi (T), ⁴% N/T. **Svenska.** ¹Lösning, ²Observerat (O), ³Förväntat (E), ⁴% O/E.

Specificity

Compound ¹	Amount Added ²	% Cross reactivity ³
Alpha2-Macroglobulin	350 mg/dL	ND
Angiotensinogen	0.3 µmol/L	ND
Angiotensin 1	100 pg/mL	ND
Angiotensin 2	100 pg/mL	ND
SHBG	170 nmol/L	ND
HCG	60,000 mIU/mL	ND
AFP	400 IU/mL	ND
Prolactin	200 ng/mL	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2015 Siemens Healthcare Diagnostics.
All rights reserved.

Made in: UK

Siemens Healthcare
Diagnostics Products Ltd.
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

0088

when used for Trisomy 21 testing
ved brug til Trisomy 21 testning
kasutatuna trisoomia 21
analüüsimiseks
lietojot 21. hromosomas trisomijas
testçðanai
kai naudojama trisomijos 21 tyrimui
når brukt til Trisomy 21-testing
når använd för trisomi 21-tester

2015-07-16

PINL2KPC – 5 {29}

Changes in this Edition:

cc#EU22593: 1) Removed references to discontinued products. 2) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS pictograms/signal words/hazard and precautionary codes and statements, offending chemicals and sources. 3) Under References, added citations from the potential biohazard statement. 4) Changed "Origin: UK" to "Made in: UK". 5) Updated the manufacturer address to include Glyn Rhonwy. 6) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

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