# IMMULITE® 2000 Multi-Diluent 1

---

## English

---

**Intended Use:** For the on-board dilution of patient samples exceeding the calibration range of the indicated IMMULITE 2000 and IMMULITE 2500 assays. See the IMMULITE 2000 and IMMULITE 2500 Operator's Manual for an explanation of on-board dilution features.

Catalog Number: **L2M1Z**

**Use with IMMULITE 2000:** Cortisol (L2KCO), DHEA-SO<sub>4</sub> (L2KDS), PAP (L2KPA), Progesterone (L2KPW), Total Testosterone (L2KTW), Total T3 (L2KT3), Total T4 (L2KT4) and Troponin I (L2KTI).

**Use with IMMULITE 2500:** Cortisol (L5KCO), DHEA-SO<sub>4</sub> (L5KDS), PAP (L5KPA), Progesterone (L5KPW), Total Testosterone (L5KTW), Total T3 (L5KT3), Total T4 (L5KT4) and Troponin I (L5KSTI).

## Materials Supplied

One vial of concentrated (ready-to-use) processed, normal human serum, containing undetectable to low levels of cortisol, DHEA-SO<sub>4</sub>, PAP, progesterone, testosterone, total T3, total T4, and troponin I, with preservative. Storage: 30 days (after opening) at 2–8°C or for 6 months (aliquotted) at –20°C. Dispose of in accordance with applicable laws.

**L2M1Z:** 25 mL.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2M1Z:** 3 labels.

## Warnings and Precautions

For *in vitro* diagnostic use.

Multi-Diluent 1 is not intended for use in linearity or sensitivity studies, but rather for the dilution of patient samples that exceed the calibration range of the assay.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were

tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

## Technical Assistance

For technical assistance, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.
©2010 Siemens Healthcare Diagnostics Inc. All rights reserved.