CV-190

Video processing powering advanced endoscopy.

Main features

- 1.3
- NBI (Narrow Band Imaging) in EVIS EXERA III 190 series scopes provides twice the viewable distance of EVIS EXERA II 180 series scopes and offers much greater contrast between blood vessels and mucosa
 - The CV-190 video processor in combination with all HQ scopes contains the necessary technology for dual focus imaging, which delivers an optimal view in two modes: near focus or normal focus
 - The newly designed waterproof one-touch connector enables a one-step connection to the light source and does not require a separate scope cable for the video processor
 - New and improved image processing delivers sophisticated image quality via enhanced colour reproduction, minimised image noise and reduced halation
 - The pre-freeze function automatically selects the clearest still image, saving time
 - Compatible with the EVIS 100/130/140/150, EVIS EXERA 160, EVIS EVERA II 180, EVIS EXERA III 190 and GI/BF/VISERA series scopes
 - 16:9 and 16:10 output for an HDTV monitor is available; compatible with analogue, HD-SDI and DVI output
 - Link connection to peripheral devices avoids complicated cable connections and accelerates transmission speed
 - OLYMPUS documentation system enhances networking expandability
 - Picture-in-picture and index function effectively enhance your observation
 - Portable memory is compatible – standard for data management; simply connect and upload
 - Supports DV output to compatible documentation devices

Specifications

Power supply	Voltage	100-240 V AC (NTSC)/220-240 V AC (PAL); within ±10%
	Frequency	50/60 Hz; within ±1 Hz
	Electric power consumption	150 VA
Size	Dimensions (W x H x D)	370 x 85 x 455 mm; 382 x 91 x 489 mm (maximum)
	Weight	10.7 kg
Classification (medical electrical equipment)	Type of protection against electric shock	Class I
	Degree of protection against electric shock of applied part	Depends on applied part. See applied part (camera head or videoscope).
	Degree of protection against explosion	The video system center should be kept away from flammable gases.
Observation	Analogue HDTV signal output	Either RGB or YPbPr (1080/50i: PAL) output can be selected. - 1.2.2
	Analogue SDTV signal output	VBS composite, Y/C and RGB (576/50i: PAL); simultaneous outputs possible. - 1.2.1
	Digital signal output	HD-SDI (SMPTE 292M), SD-SDI (SMPTE 259M), DV (IEEE 1394) and DVI (MUXGA, 1080p or SXGA) can be selected. - 1.2.3
	White balance adjustment	White balance adjustment is possible using the white balance button on the front panel.
	Standard colour chart output	The 'colour bar' or the '50% white' screen can be displayed.
	Color tone adjustment	The following colour tone adjustments are possible using the colour tone level adjustment button and colour tone selector button on the keyboard. - Red adjustment: ±8 steps - Blue adjustment: ±8 steps - Chroma adjustment: ±8 steps
	Automatic gain control (AGC)	The image can be electronically amplified when there is inadequate light due to the distal end of the endoscope being too far from the object.
	Contrast	- N (Normal): Normal image - H (High): The dark areas are darker and the bright areas are brighter than in the normal image. - L (Low): The dark areas are brighter and bright areas are darker than in the normal image.
	Iris	The auto iris modes can be selected using the 'iris mode' switch on the front panel. - Auto: The brightness is adjusted based on the brightest part of the central part and the average brightness of the periphery part. - Peak: The brightness is adjusted based on the brightest part of the endoscopic image. - Average: The brightness is adjusted based on the average brightness of the endoscopic image.
	Image enhancement setting	Fine patterns or edges in the endoscopic images can be enhanced electrically to increase the image sharpness. Either the structural enhancement or edge enhancement can be selected according to the user setup. - Structural enhancement: Enhancement of contrast of the fine patterns in the image. - Edge enhancement: Enhancement of edges of the endoscopic image.
Switching the enhancement modes	Three levels of enhancement can be selected using the image enhancement mode button on the front panel.	
Image size selection	The size of the endoscopic image can be changed using the IMAGE SIZE key on the keyboard.	
Freeze	An endoscopic image is frozen using an endoscope or the FREEZE key on the keyboard.	
Switching the method of freezing the endoscopic image	Pre-freezing: The image with the least blur is selected from the images captured in the set time period before the freeze operation and is displayed.	
Fog-free function	When a compatible endoscope is connected to the video system center, the fog free function can be used.	
Endoscope's remote switches function	The functions of the remote switches on the endoscope can be set in the user settings.	
Reset to defaults	The following settings can be reset to their defaults using the reset button on the front panel. - Colour tone - Iris mode - Image enhancement mode - Colour enhancement mode - Optical-digital observation - Image size - Contrast - Freeze - Release index - Electronic zoom - Optical-digital observation - Arrow pointer - Stopwatch - Characters on screen - PIP/POP	
Remote control	The following ancillary equipment can be controlled (specified models only). - Monitor - DVR - Video printer - Image filing system	
Documentation	Patient data	The following data can be displayed on the monitor using the keyboard. - Patient ID - Patient name - Sex - Age - Date of birth - Date of recording (time, stopwatch) - Comments
	Displaying the record state	The recording state of the following ancillary equipment can be displayed on the monitor. - Portable memory and internal buffer - DVR - Video printer - Image filing system
	Displaying the image information	The following data can be displayed on the monitor. - Structure enhancement level - Edge enhancement level - Zoom ratio - Color mode - Focus
	Advanced registration of patient data	Up to 50 patient data sets can be registered. - Patient ID - Patient name - Sex and age - Date of birth - 1.10
Portable memory	Media	MAJ-1925 (OLYMPUS)
	Recording format	- TIFF: no compression - JPEG (1/5): approx. 1/5 compression - JPEG (1/10): approx. 1/10 compression
	Number of recording images	- TIFF: approx. 227 images - JPEG (1/5): approx. 1024 images - JPEG (1/10): approx. 2048 images
Memory backup	User settings	Up to 20 user settings can be registered.
	Memorisation of selected setting	The following settings are stored even after the video system center is turned off. - Colour tone - Iris mode - Enhancement - Colour enhancement mode - Contrast - AGC - Colour mode - White balance
	Lithium battery	Life: 5 years

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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4 sk.

CLV-190

Video Xenon light source - one-touch access to HDTV and NBI.

CLV-190

Main features

- NBI (Narrow Band Imaging) in EVIS EXERA III 190 series scopes provides twice the viewable distance of EVIS EXERA II 180 series scopes and offers much greater contrast between blood vessels and mucosa thanks to an improved lamp design and signal processing
- The newly designed waterproof one-touch connector allows a one-step connection to the light source and does not require a scope cable
- Considerable reduction in operating noise thanks to a redesigned fan
- Link connection to peripheral devices avoids complicated cable connections and accelerates transmission speed
- Automatic light adjustment to achieve ideal illumination for observation with each scope

Specifications

Power supply	Voltage	100–240 V AC; within ±10%
	Frequency	50/60 Hz; within ±3 Hz
	Consumption electric power	600 VA
Size	Dimensions (W x H x D)	370 x 150 x 476 mm (standard) 390 x 162 x 551 mm (maximum)
	Weight	19 kg
	Illumination	Examination lamp
	Average lamp life	Approximately 500 hours of continuous use. (With intermittent use, the lamp life may vary slightly.)
	Ignition method	Switching regulator
	Brightness adjustment	Light-path diaphragm control
	Cooling	Forced-air cooling
	Intensity mode	Normal or high intensity
	NBI observation	Available
	Colour conversion	Possible using special-purpose filter
	Emergency lamp	Halogen lamp (within mirror) 12 V 35 W
	Average emergency lamp life	Approximately 500 hours
Automatic brightness adjustment	Automatic brightness adjustment method	Servo-diaphragm method
	Automatic exposure	17 steps
Air feeding	Pump	Diaphragm type pump
	Pressure switching	4 levels available (OFF, low, medium, high) - 2.3
Water feeding	Method	Feeds water by pressurising the detachable water container with air.
Indicators on front panel	Emergency lamp	Indicates absence of emergency lamp, disconnection and use of emergency lamp.
	NBI	When the NBI observation is enabled, the NBI indicator lights up. - 2.2
	PDD	When the PDD observation mode is enabled, the PDD indicator lights up.
Setting memory		Settings (except filter setting) are stored even when the light source is off.
Classification (medical electrical equipment)	Type of protection against electric shock	Class I
	Degree of protection against electric shock of applied part	Depends on applied part. See applied part (camera head or videoscope).
	Degree of protection against explosion	This instrument should be kept away from flammable gases.

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CF-H190L/I

A powerful tool for routine colonoscopy – standard with HDTV and RIT.

CF-H190L/I

Main features

RIT (Responsive Insertion Technology)

RIT combines three proprietary insertion tube technologies: HFT (High Force Transmission), PB (Passive Bending) and variable stiffness. These technologies are designed to facilitate complete colonoscopies by improving scope handling, insertability and ergonomics. PB helps EVIS EXERA III 190 series scopes move through acute bends in the colon. HFT provides improved operator control for both pushing and twisting manoeuvres. Variable stiffness allows the physician to adjust the rigidity of OLYMPUS scopes as needed, by simply turning an adjustment ring on the scope control section.

HDTV image quality

3.7 The CF-H190 features HDTV. Its close-up capability makes maximum use of a wide 170° angle of view and high-definition NBI images to support detailed observation. What's more, the distal end diameter is just as slim as that of the previous model.

HDTV
Image courtesy of Jean-Francois Rey, MD

NBI (Narrow Band Imaging)


NBI in EVIS EXERA III 190 series scopes provides twice the viewable distance of EVIS EXERA II 180 series scopes and offers much greater contrast between blood vessels and mucosa. The greatly improved performance of NBI opens up exciting new clinical applications and reinforces NBI's position as the standard of care for GI endoscopy.

NBI
Image courtesy of Jean-Francois Rey, MD


Waterproof One-touch Connector

A new connector design minimises the effort required for set-up prior to and in between cases. In addition, it is fully submersible and eliminates the need for a water-resistant cap and the associated risk of an expensive repair due to accidental immersion.

Specifications

Optical system	Field of view	170°	3.2
	Direction of view	Forward viewing	
	Depth of field	2-100	3.3
Insertion section	Distal end outer diameter	13.2 mm	3.16
	Distal end enlarged		

Instrument channel	Insertion tube outer diameter	12.8 mm	3.5
	Working length	L: 1680 mm, I: 1330 mm	3.8
	Channel inner diameter	3.7 mm	3.17
	Minimum visible distance	3.0 mm from distal end	
Direction from which endotherapy accessories enter and exit the endoscopic image			


Bending section	Angulation range	Up 180°	3.4.1
		Down 180°	
		Right 160°	3.4.2
		Left 160°	
Total length	L: 2005 mm I: 1655 mm		
Compatible EVIS EXERA system	Video system center OLYMPUS CV-190 Xenon light source OLYMPUS CLV-190		

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
DUAL-LOCKING

V
.SYSTEM™

Unique dual guidewire locking mechanism for fast, secure short guidewire exchange, enhancing flexibility and reliability


Setting a New Standard for Ease of Use, Security, and Flexibility, the Unique Dual Locking Mechanism Supports Enhanced ERCP Efficiency

4.8

Dual Locking Mechanism Securely Locks 0.025" and 0.035" Guidewires

Completely redesigned to ensure greater reliability and flexibility, the TJF-Q180V's dual locking mechanism is optimized to exploit the reactive force of the guidewire. The forceps elevator has been modified to broaden the range of scope positions in which the Side Lock mechanism.

4.8.1

Center Lock

With its firmer grip and a design that takes advantage of the guidewire's reactive force, the Center Lock section, the V-Groove, on the forceps elevator has a reconfigured groove shape that locks the guidewire more securely than ever.

When the papilla is in the green area* as shown in the figure on the right, the Center Lock will be used.

Side Lock

4.8.2

A new Side Lock section has been added to the side of the forceps elevator to increase guidewire locking flexibility. It physically fixes the guidewire even when it is positioned to the side of the forceps elevator.

When the papilla is in the blue area* as shown in the figure on the right, the Side Lock will be used.


*Illustrations are colored for greater clarity. The monitor screen is not actually colored like these figures.


Excellent Image Quality

The TJF-Q180V incorporates a high-resolution CCD that delivers the sharp, clear images Olympus is known for, displaying them in a large screen size for easy viewing. Advanced Narrow Band Imaging™ capability is also supported with this scope.

4.1

Normal image

NBI image

Conventional elevator wire channel port

New sealed elevator wire channel port

Easier Cleaning (Unique Fixed Distal End Design)

In addition to the clinical performance benefits of the TJF-Q180V's locking mechanism, the elevator wire channel port is now sealed so separate cleaning is no longer necessary. The result is faster, easier cleaning that makes scope reprocessing more efficient.

Main Features


4.1 {

- Unique dual locking mechanism for efficient short guidewire exchange in combination with dedicated V-System™ ERCP devices.
- Superb image quality with exceptionally sharp, crisp details, and large display size.
- Narrow Band Imaging™ enhances observation of the mucosa and capillaries.
- Ergonomically designed grip enhances scope maneuverability, while easy-to-access controls and user programmable switches improve operability.
- Slim 11.3 mm insertion tube.
- Wide 4.2 mm diameter channel.
- 4-way angulation (120° up, 90° down, 110° right, and 90° left) facilitates approach to the papilla of Vater.
- Fully compatible with the CV-160 and 140; use of CV-180 enables the full feature set.
- Scope ID function stores individual scope information in the built-in memory chip and displays it on the monitor, facilitating endoscopy suite management.

Specifications

Optical System	Field of view	100°	- 4.2
	Direction of view	Backward side viewing 5°	
	Depth of field	5 to 60 mm	
Distal End	Outer diameter	13.7 mm	- 4.5
	Outer diameter	11.3 mm	
Insertion Tube	Outer diameter	11.3 mm	- 4.4
Bending Section	Range of distal end bending	Up 120°, Down 90°, Right 110°, Left 90°	- 4.3
Working Length		1,240 mm	- 4.7
Total Length		1,550 mm	
Instrument Channel	Inner diameter	4.2 mm	- 4.6
	Minimum visible distance	10 mm	
	Endotherapy accessory entrance/exit position in field of view		

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OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of sales and distribution locations visit:
www.olympus.com

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
EVIS EXERA III

GIF-H185

Gastrointestinal videoscope – outstanding HDTV imaging for routine endoscopy.

Main features

HDTV image quality

With the new EVIS EXERA III system, HDTV image quality delivers high-definition observation capabilities, even with the slim scope design.

NBI (Narrow Band Imaging)

NBI in EVIS EXERA III 185 series scopes provides twice the viewable distance of EVIS EXERA II 180 series scopes and offers much greater contrast between blood vessels and mucosa. The greatly improved performance of NBI opens up exciting new clinical applications and reinforces NBI's position as the standard of care for GI endoscopy.

Close focus

Close focus enables you to obtain an enlarged, close-up image simply by moving the scope tip as close as 2 mm from the mucosa.

Slim design

This scope offers an excellent balance between size and performance, with HDTV image quality in a slim 9.2 mm diameter size.

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
Specifications

Optical system	Field of view	140° - 512
	Direction of view	Forward viewing
	Depth of field	2-100 mm - 513
Insertion section	Distal end outer diameter	9.2 mm - 515
	Distal end enlarged	
Instrument channel	Insertion tube outer diameter	9.2 mm - 516
	Working length	1030 mm - 518
	Channel inner diameter	2.8 mm - 514
	Minimum visible distance	3.0 mm from the distal end
	Direction from which endotherapy accessories enter and exit the endoscopic image	

Bending section	Angulation range	Up 210°
		Down 90°
		Right 100°
		Left 100°
Total length	1350 mm	
Compatible EVIS EXERA system	Video system center OLYMPUS CV-190	
	Xenon light source OLYMPUS CLV-190	

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




HIGH DEFINITION LCD MONITOR

OEV262H

26-Inch LCD Monitor having the advanced display technology to realize the full potential of the Olympus Endoscope System

HIGH DEFINITION LCD MONITOR OEV262H

Main Features

6.1

6.2

6.5

- 26-inch Full HD LCD Panel with a 16:9 aspect ratio, high brightness, high contrast, and high gradation images.
- FLIP function (such as mirror and 180° rotation) which provide suitable display pattern and monitor layout for the procedure.
- Multi-modality display capability, including Picture-In-Picture (PIP), Picture-Out-Picture (POP), and Clone Out^{*1} with various image size combinations.
- Various Input/ Output terminals, including 3G/HD/SD SDI (x2), DVI-I (x2), HD15, Y/C, and Composite.
- Eco-friendly specifications by low power consumption, various power saving mode, lightweight and thin body compared to predecessor^{*2}.
- A.I.M.E (Advanced Image Multiple Enhancer)^{*3} which supports to enhance the image quality of the OLYMPUS endoscope system.

6.5.1

6.5.2

*1 Clone Out is the function that outputs the HD video signal (1080i/1080p) as displayed including PIP/POP image to the second monitor and record device.
 *2 HIGH DEFINITION LCD MONITOR OEV261H.
 *3 A.I.M.E. (Advanced Image Multiple Enhancer) is the function that improves the resolution with saving the noise gain up.

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














MOBILE WORKSTATION
WM-NP2

Optional Accessories

Customisation of the workstation is important when considering the application for which it will be used - the WM-NP2 is no exception. A wide range of accessories is available as follows:

MAJ-1638 Scope Pole Kit *		
MAJ-1639 CO ₂ Cylinder Holder (140 mm)		
MAJ-1640 Sliding Keyboard Tray *		
MAJ-1641 LCD Monitor Arm *		
MAJ-1642 IV Pole		
MAJ-1643 Drawer Unit *		
MAJ-1650 CO ₂ Cylinder Holder (205 mm)		
MAJ-1652 CO ₂ Cylinder Holder (Double) *		
MAJ-1654 Equipotential Strip		
MAJ-1655 OEV 26" Monitor Power Supply Holder		
MAJ-1656 Footswitch Holder		
MAJ-1657 Keyboard Arm - Side Mounted		
MAJ-1658 Universal Storage Container		
MAJ-1659 Equipotential Support Rail		
MAJ-1660 Cable Retaining Hook		
MAJ-1662 Nurses Control Arm		
MAJ-1663 EUS Arm Mount Kit		
MAJ-1664 Nurses Control Panel Arm - Long		
MAJ-1666 Side Shelf		
MAJ-161 Rail Clamp		

* provided as standard in some kits

Specifications

WM-NP2

Dimensions: 1826 mm (H) x 700 mm (D) x 650 mm (W)

Height to Top Tray work surface 1169 mm

Height with 26" LCD monitor in highest position 2073 mm

Weight: Maximum 85 kg (unladen) including fitted transformer

Voltage

The WM-NP2 is supplied with a separating transformer based on one of the following voltage options

MAJ-1646 WM-T2 Transformer 100 V

MAJ-1647 WM-T2 Transformer 110-120 V

MAJ-1648 WM-T2 Transformer 220-240 V


Construction

Workstation

Hoop: Extruded profiled aluminium

Shelves: Mild steel

Top tray: Injection moulded ABS

MOBILE WORKSTATION
WM-NP2
New Workstation for EVIS EXERA III

Specifications, design and construction are subject to change without any notice or obligation on the part of the manufacturer.

For a complete listing of sales and distribution locations visit www.olympus.com

Issue 1, May 2012

Printed in UK. K1 5036016/512

The Second Generation WM-NP2 Series Mobile Workstation has been designed specifically for the EVIS EXERA III Video System, incorporating an innovative and practical design.

Sleek and sophisticated in design, practical in use

Olympus' R&D Department has 40 years' experience of workstation design and has developed practical platforms to support all areas of surgical and medical endoscopy, culminating in this the most advanced workstation yet, the WM-NP2.

Endoscopic procedures have continued to develop and the interface between user equipment remains critical to provide a safe and time efficient procedure. As a result the importance of the workstation as the 'imaging' and equipment platform is crucial in providing access to all devices and particularly the viewing monitor.

With the introduction of a gas assisted LCD Monitor Arm, the ability to manipulate the monitor for the optimum viewing angle improves functionality and reduces strain on the clinician undertaking the procedure.

The WM-NP2, setting a new standard in endoscopic workstation design.

Storage

The WM-NP2 will accommodate a wide range of electro-medical devices to meet the requirements of either routine or advanced endoscopic procedures.

Monitor position

The all new LCD Monitor Arm allows the operator to optimise the viewing position with ease, due to the integral gas spring and tensioned joints providing excellent balance in use.

Drawer pack

The drawer pack facilitates safe storage of a variety of devices, providing a practical addition to the workstation. Multiple drawer packs can be fitted to a single workstation.

Cable management

Optimised cable management is achieved through a range of cable channels integral to the workstation hoop with the addition of dedicated cable winders at the rear of the workstation shelves.

Electrical safety

An Olympus separation transformer is fitted as standard providing the required level of electrical safety. Output from the separation transformer is controlled using the central on/off switch, allowing all equipment on the workstation to be 'powered up' simultaneously.

Keyboard tray


Depending upon set configuration, the sliding keyboard tray is simple to use and enables easy access for entering demographic data.

Manoeuvrability

With increasing volumes of equipment being stored on the workstation, ease of movement is critical. The WM-NP2 has handles incorporated into the top tray, on the front, side and rear facing surfaces.

'Clip in' Accessories

A completely new concept, a range of 'clip in' accessories are provided, from water bottle clips to cable holders allowing users to customise their workstation.

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SUCTION PUMP

KV-6

The KV-6 suction pump is the successor to the KV-5 pump. It is designed for endoscopic aspiration and for non-endoscopy applications. The unique flow change buttons allow a fast vacuum build-up with a single touch. It has a "whisper mode" for quiet operation and "turbo mode" for faster suction. The pump connects to an extensive fluid collection system.

Key Benefits

- **High Suction Performance**

The KV-6 suction pump has vacuum capabilities of up to -95 kPa. There are three pre-set speed settings (40 l/min, 50 l/min, and 60 l/min), which enables the operator the flexibility to control the rate of vacuum build-up at the touch of a button.

- **Cleaning and Hygiene**

- Single-use collection system allows easy and safe disposal of clinical waste.
- Touch-sensitive buttons are easy to clean.

- **Compatibility**

Compatible with existing KV-5 accessories and Olympus workstations.

The KV-6 suction pump is ideal for use either mounted on a procedure trolley or as a stand-alone unit. The variable suction control allows adjustment to meet individual requirements. A high-performance pump and motor combination ensure that the unit is "whisper quiet" for any working environment.

The KV-6 may be adapted to suit individual requirements by the use of optional components, to allow a single-use or reusable suction system. The single-use liners are conveniently supplied with solidifier for safe and simple use and disposal.

Product Details

Compact Design for Convenient Use

The KV-6's compact and easy-to-use design means that it is suitable for use in office space or the operating theater.

Simple and Easy to Operate with Clearly Labelled Controls

A clear control panel ensures that the KV-6 suction pump can be used safely and efficiently by all staff.

Trolley-Mountable for Convenience

The KV-6 is ideal for use either mounted on a procedure trolley or as a stand-alone unit for maximum flexibility of use.

82

Use with 1.5- or 2.5-Liter Jars and Liners

The KV-6 can be used with either 1.5- or 2.5-liter jars, enabling its use in a variety of procedures.

Suction Jar Mountable on Olympus Workstations


As the suction jar can be mounted on Olympus workstations, it is compatible with the Olympus range of devices and equipment for the office or operating theater.

Overview of Accessories

Suction Tubing

K10026130
with 2 coupling pieces (1 x 90°),
7 x 12 mm,
60 cm

K10026804
with 2 coupling pieces,
7 x 12 mm,
60 cm

Patient Connecting Tubing

K7503432
50 pack

Disposable Liners

K10021894
K10021895
1.5 liter

K10021897
K10021898
2.5 liter

Reusable Lids

K10026127
0.25 liter

K7503760
1 & 2 liter

Jars for Liners

K10021893
1.5 liter

K10021896
2.5 liter

Reusable Jars

K10026128
0.25 liter

K7503762
1 liter


K7503763
2 liter

Filters

K10027321
Disposable bacteria/odor filter
10 pack

K10027322
Disposable bacteria filter
10 pack

K10026130
K10026127

K10026128

Clamp Holder

K10002819
MAJ-191

K10021793
MAJ-1659
K10001115
MAJ-190

K7503556
WM-60 / D-60

Two or More Jars

Disposable collection system

K10021885
MAJ-1673

K10021188
10 Pack

K10021878
MAJ-1672

Disposable/reusable collection system

K10021882
MAJ-1674

Specifications

Classification (Electromedical Equipment)	Standards Compliance	Complies with EN ISO 10079-1 and EN IEC 60601-1 and UL 60601-1.		
	Electromagnetic Compatibility	This product complies with the requirements of EN IEC 60601-1-2 for emissions and immunity, and as such its operation is unlikely to be affected by, or cause interference with, equipment meeting appropriate EMC standards. As a precaution, equipment which may be sensitive to interference outside the limits specified by EN IEC 60601-1-2 should not be placed in close proximity to the KV-6.		
	Type of Protection against Electrical Shock	Classification according to EN IEC 60601-1 and UL 60601-1: Class 1 with Type CF applied part.		
	Degree of Protection against Electrical Shock	In accordance with EN IEC 60601-1 and UL 60601-1, the KV-6 is marked with the symbol to indicate the provision of an adequate degree of protection against electrical shock and that it has an applied part isolated from all other parts of the equipment.		
	Degree of Protection against Explosion	None: the Olympus KV-6 must NOT be used in the zone of risk of flammable anesthetic gases.		
	Mode of Operation	Continuous.		
Regulatory Status	European Economic Area (EEA)	CE: This mark on the product indicates compliance with Directive 93/42/EEC relating to medical devices, Class IIa.		
Dimensions and Weight	Dimensions	Height: 210 mm	Width: 305 mm	Depth: 375 mm
	Weight	9.3 kg		
Power Switch	Marking	The power switch is marked. A white LED in the center signalizes it is turned on.		
Power Requirements	Power Supply	100–240 V~	Fusing	2 x T1.6 AH
	Frequency	50/60 Hz	Power Rating	120 W
	Marking: The mark ~ on the product indicates the requirement for an AC power supply.			
Environmental Conditions	Ambient Temperature	Operational:	+5°C to +40°C (+41°F to +104°F)	
		Storage:	-20°C to +50°C (-4°F to +122°F)	
	Relative Humidity	Operational:	30–75%	
		Storage:	20–95%	
Atmospheric Pressure	Operational/Storage:	70.0 to 106 kPa		
Fluid Ingress	In accordance with EN IEC 60601-1 and UL 60601-1, the product is marked IP21 to indicate that it is provided with an enclosure that prevents entry of such an amount of falling liquid as might interfere with its safe and satisfactory operation when correctly positioned.			
Resistance to Chemicals	The external surfaces of the Olympus KV-6 are resistant to: 2% aqueous neutral detergent, 70% ethyl alcohol, isopropyl alcohol, water.			
Pump Specification	Nominal Vacuum:	95 kPa	Nominal Free Air Flow Rate:	50 l/min, adjustable to 40 or 60 l/min.
	In accordance with ISO 10079-1, the unit is marked with the words "HIGH VACUUM, HIGH FLOW" to indicate it has a vacuum of at least 60 kPa in 10 seconds and a free air flow rate of 20 l/min or greater.			
Thermal Cut-Out	The temperature is observed by a sensor in the housing. If the temperature in the housing exceeds 80°C (176°F), the software stops the motor and turns on the yellow LED. It is possible to restart the motor as soon as the temperature in the housing drops below 70°C (158°F).			
Single-Use Suction Liner	To prevent fluid from being drawn up to the pump the single-use suction liner contains a mechanism to shut off the flow when the liner is overfilled. This mechanism is activated when contact is made with aspirated fluid.			
Suction Jar	The suction jars can be autoclaved up to 135°C or disinfected. The jar is impact-resistant.			
Suction Jar Lid	The suction jar lid clips and where applicable the float mechanism are autoclavable up to 135°C. The inlet connection (from the endoscope) is marked with the word PATIENT, and the outlet connection (to the pump) is marked with the word VACUUM.			
Tubing	Use only tubing specified in the KV-6 Instructions for Use, otherwise equipment malfunction may result.			
Vacuum Gauge	Range:	0 to 100 kPa	Accuracy:	±2.5%

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

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OLYMPUS FLUSHING PUMP OFF-2

9.1


With its powerful but controllable peristaltic pump the OFP-2 will irrigate fluid down either the instrumentation or auxiliary water channel of Olympus Gastrointestinal EVIS ACTERA, EXERA, LUCERA, SMARTAGE and EUS endoscopes.

The OFP-2 facilitates two main functions:

- The washing of gastric and colonic mucosa resulting in improved visualisation, diagnosis and therapy during endoscopic procedures
- To aid the use of trans-endoscopic ultrasound probes by providing rapid filling of the organ to be examined

Main Unit:

- Designed to prevent over pressurisation of Olympus endoscopes
- Microprocessor control ensures accurate and precise, flow rate delivery
- Ergonomically designed touch buttons for simple use
- Bright LED display for easy verification of selected settings


Practicality:

- Use of the MAJ-1606 irrigation adaptor allows either irrigation or a combination of irrigation and the passage of an EndoTherapy device through the endoscopes instrumentation channel

Safety Features:

- Pump Head automatically cuts out if opened during operation preventing accidental user injury
- Standby mode allows quick and safe exchange of water channel tubes without the need to operate the mains switch
- A 20 second automatic cut-off timer allows controlled irrigation



9.4

Water Container:

- Holds two litres with markings to indicate fluid levels
- A specially designed lid guides the water tubing to optimise fluid usage
- The water container and lid are autoclavable for effective sterilisation
- Washer disinfector compatible

9.8.5.7

Size and functionality:

- Improved functionality allows use with either the instrumentation or auxiliary water channel, including EUS ultrasound endoscopes for the first time
- Compact design for essential space saving on workstations
- Operates via either remote controlled* scope operation or foot switch

9.5


9.6

* requires the use of MAJ-920 remote control cable

Consumables:

- The MAJ-1607 and MAJ-1608 channel water tubes and the MAJ-1606 adaptor are supplied sterile for immediate use
- The MAJ-1608 auxiliary channel water tube is designed for one single days use
- Instrument channel adaptor and MAJ-1607 water tube are single patient use items

MAJ-1606


Single patient use

Supplied sterile
Instrument channel adaptor

Pack 100: Part number - K10007072
Pack 10: Part number - K10016091

MAJ-1607

Single patient use

Supplied sterile
Instrument channel water tube

Pack 50: Part number - K10001146
Pack 10: Part number - K10016136

MAJ-1608

Single days use

Supplied sterile
Auxiliary channel water tube

Pack 50: Part number - K10001147
Pack 10: Part number - K10016135

9,8.52

Accessories:

MAJ-1603

Pack of three

Water container two litre
Part number - K10007071

MAJ-855

Reusable

Supplied non sterile
Auxiliary water tube
Part number - MAJ-855

MAJ-920

Remote control cable
Part number - MAJ-920

Specification:

Dimensions:	200(W) x 173(H) x 385(D) mm	Weight:	4 kg (with empty water container)
Maximum Flow Rate:	Instrument channel: 750 ml/min - 9,3 Auxiliary water channel: 230 ml/min - 9,2		
Safety:	Complies with EN/UL/IEC60601-1	Regulatory:	CE marked and ETL Listed
Power Requirements:	Power Supply: 100-240 V - Fusing: 2 x T2.5 AH 250 V	Frequency:	50/60 Hz,
		Power rating:	100 VA

Olympus reserves the right to alter the specification without notice.

OLYMPUS KEYMED GROUP
 KeyMed House, Stock Road, Southend-on-Sea,
 Essex SS5 5QH, United Kingdom
 Telephone: +44 (0)1702 616331 Fax: +44 (0)1702 466577
 e-mail: info@olympus.co.uk www.olympus.co.uk

OLYMPUS MEDICAL SYSTEMS CORP.
 Shinjuku Mainbldg. 3-1, Nishi-Shinjuku 2-chrome, Shinjuku-ku, Tokyo 163-0914,
OLYMPUS EUROPA HOLDING GMBH
 Prommes/Guedh delving Wandestraße 14-16, 20097 Hamburg, Germany
 (Letters Postfach 10 49 06, D-20034 Hamburg, Germany)
OLYMPUS AMERICA INC.
 3500 Corporate Parkway, PO Box 610, Center Valley, PA 19934-0610, USA
OLYMPUS SINGAPORE PTE LTD.
 451B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 245573, Singapore

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