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*(Devices in class Ila, I Ib or III) according to EC Certificate
Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)
issued by Notified Body TÜV SÜD Product Service GmbH, Ridlerstrasse 65,
80339 Munich, Germany, CE 0123*

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to which this declaration relates, comply with all relevant provisions of Council Directive 93/42/EEC. The products comply with the standards and/or other normative documents as listed in the applicable technical documentation.

Ulm, den 09.10.2013

A handwritten signature in black ink, appearing to read "R. Bohnacker".

Rolf Bohnacker
Representative Quality Management

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EN Instructions for use

Patient tubing for ulrich medical CT/MRI injectors

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Highlighting and symbols

Highlighting	Meaning
	Hazards to persons. Non-compliance can cause death or serious injuries to persons and damage to the product.
→	Handling instructions.

Symbols on the product and packaging

Icon	Description	Icon	Description
	Caution		Sterilized using ethylene oxide

Manufacturer's CE mark

Do not use if package is damaged

REF	Product number	Do not use tubing if the protective caps are not properly attached
	Lot code	Pyrogen-free
	Use by date	Keep dry
	Manufacturer	Temperature limit
	Consult the instructions for use	Relative humidity limit
	Do not re-use	Protect from sunlight

Intended use

The patient tubing for ulrich medical CT/MRI injectors is used for venous-side administration of contrast agents and physiological saline solutions (NAC) during a CT/MRI examination. The patient tubing creates the connection from the ulrich medical pump tubing (device side) to the patient.

Product description

The patient tubing consists of tubing with two protective caps and has a male Luer lock connector with integrated check valve on the patient side. The pump tubing is connected to the female Luer lock connector. A second check valve is integrated in the tubing. The check valves are self-closing. They only open when a medium flows in the defined direction. If back pressure is exerted or if there is no media flow, the valves remain closed.

A	Device side with pump tubing
B	Patient side
1	Female Luer lock
2	Protective caps for tubing ends
3	Male Luer lock with first check valve
4	Second check valve
5	Media flow direction

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Art. no.	Designation of the patient tubing
XD 2030	Patient hose 1.20 m
XD 2035	Patient hose 1.50 m
XD 2040	Patient hose 2.50 m
XD 2045	Patient hose 3.20 m

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- Do not reuse product.
- Do not resterilize product.
- Do not use product from open or damaged sterile packaging.
- Discontinue using product if the use-by date has expired.
- Do not use product that was not stored under the recommended storage conditions.

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- Before using the product, check its operability and proper condition.
- To avoid damage due to improper use, use product according to these instructions for use and the instructions for use of the injector
- Product and accessories should only be used by persons who have the required knowledge, training or experience.

Risk of injury!

The use of tubing combinations or systems that have not been approved by ulrich medical jeopardizes the safety of the patient and/or the user as well as the proper functioning of the injector.

- Only use original products approved by ulrich medical.

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Tubing extensions jeopardize the safety of the patient and risk contaminating the tubing system.

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Risk of injury due to infections and contaminations!

Improper handling jeopardizes the safety of the patient.

- Handle all sterile components carefully in order to ensure sterility.
- Do not use product if the protective caps are missing or are not attached securely.
- Replace patient tubing for each new patient.
- If new patient tubing is not connected immediately after disconnecting the patient tubing from the pump tubing: Place a sterile protective cap on the pump tubing.
- Dispose of patient tubing upon disconnecting it from the patient in order to prevent contamination of the pump tubing.

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Connect product to pump tubing

- Connect product only as directed in the instructions for use of the injector.

Performance testing

- Conduct performance testing only as directed in the instructions for use of the injector.

Operation

- Operate product only as directed in the instructions for use of the injector.

Validated processing method

Single-use product

! WARNING! Patient/user infection due to re-use and/or processing.

If re-used and/or processed, there is a risk of contamination and material fatigue. Soiling, contamination, impaired function may lead to injury, illness or death.
→ Do not process product.

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Ambient conditions

	Use	Storage
Temperature	+5°C ... +40°C	+15°C ... +25°C
Relative humidity	10% ... 90%, non-condensing	30% ... 60%, non-condensing


Disposal

- Dispose of product after use according to local regulations.

Ulrich GmbH & Co. KG (Kundendienst) Tel.: +49 (0)731 9654-0
Buchbrunnenweg 12, 89081 Ulm, Germany
E-Mail: injectors@ulrichmedical.com

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Technical specifications – tubing system for 3T™ tennessee™, mississippi™, missouri™, ohio M™ and ohio tandem™

Consumption		Length(s)	Specials
Pump tubing	Patient tubing		
1 tubing per working day (max. 24 h)	1 tubing per patient	120cm	2 check valves
		150cm	Proven resistance to pressure
		250cm	Proven CA compatibility
		320cm	Non-pyrogenic
			Latex-free
			DEHP-free
			Proven resistance to pressure
			Proven CA compatibility
			Non-pyrogenic
			Latex-free
			DEHP-free
			Particle filter
			Integrated pressure control unit
			3 accesses for reservoir containers, each with air filter

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Technical specifications – tubing system for 3T tennessee™, mississippi™, missouri™, ohio M™ and ohio tandem™

Consumption		Length(s)	Specials
Pump tubing	Patient tubing		
1 tubing per working day (max. 24 h)	1 tubing per patient	120cm	2 check valves
		150cm	Proven resistance to pressure
		250cm	Proven CA compatibility
		320cm	Non-pyrogenic
			Latex-free
			DEHP-free
			3 accesses for reservoir containers, each with air filter
			Particle filter
			Integrated pressure control unit
			Proven resistance to pressure
			Proven CA compatibility
			Non-pyrogenic
			Latex-free
			DEHP-free

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
EN Instructions for use

Pump tubing for ulrich medical CT/MRI injectors

About This Document

- These instructions for use are part of the product and they describe its safe use as intended.
- Read the instructions for use before using the product.
 - Keep the instructions for use with the product so that they are accessible at all times.

Highlighting and Symbols

Highlighting	Meaning
WARNING	Hazards to persons. Non-compliance can cause death or serious injuries to persons and damage to the product.
→	Handling instructions.

Symbols on the Product and Packaging

Icon	Description	Icon	Description
	Caution		Do not use if package is damaged
	Manufacturer's CE mark		Do not use tubing if protective cap not properly attached
	Product number		Pyrogen-free
	Lot code		Keep dry
	Use by date		Temperature limit
	Manufacturer		Relative humidity limit
	Consult the instructions for use		Protect from sunlight
	Sterilized using ethylene oxide		Do not re-use! Pump tubing intended for single use only for 24 hours.

Intended use

The pump tubing for ulrich medical CT/MRI injectors is used for venous-side administration of contrast agents and physiological saline solutions (NaCl) during a CT/MRI examination. The pump tubing creates the connection from the injector (device side) to the patient tubing. The pump tubing can be used for any number of injections within 24 hours.

Product description

On the patient side, the pump tubing consists of a tube with a protective cap which has a male luer lock connector (XD8003: with integrated check valve. The check valve is self-closing. It only opens when a medium flows in the defined direction. If back pressure is exerted or if there is no media flow, the valve remains closed). The spikes contain an integrated bacteria filter in the supply air duct.

A particle filter in the pump tubing retains foreign bodies. The pressure sensor unit controls the flow of injection volume in the pump tubing. The hygienic safety has been proven by an independent testing facility.

- A Device side with pump tubing
B Patient side with patient tubing
1 Spike for media supply for contrast agent
2 Spike for media supply for NaCl
3 Spike for media supply for contrast agent

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- 5 Tubing section for roll pump
- 6 Pressure sensor unit
- 7 Particle filter
- 8 Male Luer lock, with check valve (XD 8003), without check valve (XD 2020)
- 9 Protective caps for tubing ends
- 10 Media flow direction

Delivery contents and scope of application

Art. no.	Designation
XD 8003	Pump tubing for CT motion Injektor
XD 2020	Pump hose for XD200x Injectors

Product for single use

The product is EO sterilized and delivered in sterile packaging.

- Do not reuse product.
- Use product for a maximum of 24 h
- Do not sterilize product.
- Do not use product from open or damaged sterile packaging.
- Discontinue using product if the use-by date has expired.
- Do not use product that was not stored under the recommended storage conditions.

Safe handling

- Before using the product, check its operability and proper condition.
- To avoid damage due to improper use, use product according to these instructions for use and the instructions for use of the Injektor.
- Product and accessories should only be used by persons who have the required knowledge, training or experience.

Risk of injury!

The use of tubing combinations or systems that have not been approved by Ulrich medical jeopardizes the safety of the patient and/or the user as well as the proper functioning of the Injektor.

- Only use original products approved by Ulrich medical.

Risk of injury! Air injection, overdose and flow reduction due to pressure regulation if tubing extensions are used.

- Do not use tubing extensions (e.g. "Heidelberg extension").

Risk of injury! Air injection and spraying of contrast agent. Risk of damage to property.

If screwed in too tightly, the thread of the Luer lock connector may break and lead to open cracks.

- Close the Luer lock connector properly.

Risk of injury due to infections and contamination!
Improper handling jeopardizes the safety of the patient.

- Handle all sterile components carefully (see chart, pos. 1, 2, 3, 8) to ensure sterility.
- Do not use product if the protective caps are missing or are not attached securely.
- If new patient tubing is not connected immediately after disconnecting the patient tubing from the pump tubing: Place a sterile protective cap on the pump tubing
- Replace and dispose of the pump tubing after every working day at the latest or after 24 hours at the latest.

Risk of injury due to mixing of contrast agents with the Tandem software option!

- If a type of contrast agent other than the type previously used is loaded onto the Injektor: Change pump tubing.

Use

Preparation

- Remove the product from the packaging.
- Check product before use for damaged parts.
- Do not use product with damaged parts.

Connect product to Injektor

- Connect product only as directed in the instructions for use of the Injektor.

Performance testing

- Conduct performance testing only as directed in the instructions for use of the Injektor.

Operation

- Operate product only as directed in the instructions for use of the Injektor.
- Replace the pump tubing after every working day at the latest or after 24 hours at the latest.

Validated processing method

Single-use product

WARNUNG! Patient/user infection due to re-use and/or processing.

If re-used and/or processed, there is a risk of contamination and material fatigue. Soiling, contamination, impaired function may lead to injury, illness or death.

- Do not process product.

Technical service

Ulrich GmbH & Co. KG
Injektor Systems Customer Service
Buchbrunnweg 12
89081 Ulm / Germany
Telephone: +49 (0)731 9654-0
Fax: +49 (0)731 9654-2805
Email: injectors@ulrichmedical.com

Technical Data

Parameter	Value
Pressure tightness	20 bar max
Material	PUR, PVC (DEHP-free, latex-free)
Max. permissible duty	24 hours for multiple patients
Filling volume (XD 2020)	173 mL
Filling volume (XD 8003)	195 mL

Storage

- Store product in sterile packaging in a dry, dark and temperature-controlled room.
- Protect product in sterile packaging from environmental influences such as soil, dust, moisture, heat and sunlight.

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© 2024 ulrich medical