

Application of Titanium Clips

A56790A Clip applicator,
10 x 330 mm,
for clips medium/large A5635

- Curved
- Rotatable jaws
- Can be disassembled for reprocessing.

A5635 Clip, titanium,
medium/large,
18 x 6 pcs.


- 1 p02
- 3 p02

Use

Remove clip from the magazine.



Apply clip.



Fascial Closure

A51000A Needle,
for fascial closure

250 mm



WA51203A Needle,
for fascial closure

140 mm



WA51205A Spare cannula,
for WA51203A



Compatible Components

