

Product-Information-Sheet

Darwin's rotatable clips applicator

Easy-to-Use Reusable Rotatable Clip-Fixing Devices



- Fast and easy clip system
- For mechanical hemostasis of upper and lower GI bleeding
- Additional uses include closure of mucosal defects, marking, and anchoring
- Rotatable



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	brand_name	Area Compatibility
HX-110UR	N1085430	1	2.8 mm	2300 mm	EZ Clip	

1.1

Darwin's clips

Product-Information-Sheet

Reusable Mechanical Lithotripters ? Outer Sheaths



- Unique double sheath to aid cannulation
- Quick and effective crushing
- Injection of contrast medium (MAJ-243 only)
- System consists of a metal coil sheath, a tube sheath, a basket and a BML handle (MAJ-440)
- Compatible with emergency lithotripter (BML-110A-1)



Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	For Use With	Area Compatibility
MAJ-242 2.3	BML-3Q-1	026522	1 metal-coil sheath	4.2 mm	1950 mm	MAJ-243, MAJ-244, and MAJ-440	1

Reusable Mechanical Lithotripters ? Outer Sheaths




- Unique double sheath to aid cannulation
- Quick and effective crushing
- Injection of contrast medium (MAJ-243 only)
- System consists of a metal coil sheath, a tube sheath, a basket and a BML handle (MAJ-440)
- Compatible with emergency lithotripter (BML-110A-1)



Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	For Use With	Area Compatibility
MAJ-243 2.2	BML-3Q-1	026523	1 tube sheath	4.2 mm	1950 mm	MAJ-242, MAJ-244, and MAJ-440	1

Reusable Mechanical Lithotripter ? Baskets



- Quick and effective crushing
- Available in different basket specifications
- System consists of a metal-coil sheath, a tube sheath, a basket, and a BML handle (MAJ-440)
- Compatible with emergency lithotripter (BML-110A-1)



Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	Basket ø	For Use With	Area Compatibility
MAJ-244 2.1	BML-3Q-1	026524	1	4.2 mm	1950 mm	30 mm	MAJ-242, MAJ-243, and MAJ-440	1

Product-Information-Sheet

Reusable Eight-Wire Stone-Extraction Basket



- Unique eight-wire design for small stones and sludge
- Exceptional grasping capability ensures stone retrieval
- Bullet-shaped distal tip facilitates insertion into the bile duct
- Injection port
- Integral handle
- Compatible with emergency lithotripter (BML-110A-1)



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Basket ø (A)	brand_name	Area Compatibility
FG-301Q	026767	1	2.8 mm	1950 mm	20 mm	FlowerBasket	①

4.2 Daughardine's asthonia vialy alimony ishaulino kapselis Aerobolavojam

Reusable Stone-Extraction Baskets



- Exceptional grasping capability ensures stone retrieval
- Bullet-shaped distal tip facilitates insertion into the bile duct
- Injection port
- FG-22Q-1, -23Q-1: Compatible with emergency lithotripter (BML-110A-1), integral handle
- FG-18Q-1: To use with handle MA-479




Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Features	Basket ø (A)	Basket Wire	Area Compatibility
FG-22Q-1	026760	1	2.8 mm	1950 mm	For easy opening in narrow lumen	22 mm	Hard	①

4.1 End. Kanalo aliamubas Darbinis ilgis Kapselis astidarymo aliandubas Aerobolavojam

Reusable Triple-Lumen Sphincterotomes



- Separate lumina for guidewire, cutting-wire, and injection of contrast medium
- Stable orientation of the cutting-wire
- Distal marking system to aid sphincterotomy
- Excellent visibility under fluoroscopy
- Precurved



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Distal Tip ø	Cutting Wire (B)	Distal Tip Length (A)	Area Compatibility
KD-301Q-0725	025863	2	2.8 mm	1950 mm	4.5 Fr.	25 mm	7 mm	①

3.1 Endoskopu varakui Darbinis ilgis Distalinio gele aliamubas Dujimo vialos ilgis Distalinio gele ilgis 1/1

Product-Information-Sheet

Doughanline's biopsies enyplis

Reusable Oval-Cup Biopsy Forceps with Needle




- Needle for improved anchorage
- Elongated cups for deeper biopsies
- Fenestrated cups for larger biopsies



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Area Compatibility
FB-24K-1	026631	1	2.8 mm	1550 mm	1 1
FB-24U-1	026633	1	2.8 mm	2300 mm	1

5.4
Doughanline's griebiancias enyplis kisakio formas
Darbinis ilgis

Reusable Grasping Forceps ? Three-Prong



- For polyp and foreign-body retrieval
- Ideal for larger polyps
- Plastic-coated sheath



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Opening Width ø	Area Compatibility
FG-45U-1	026795	1	2.8 mm	2300 mm	20 mm	1 1

5.2
Endoskopo kanalas
Darbinis ilgis
Asidarymo plotis

Reusable Rat Tooth and Alligator-Jaw Rotatable Grasping Forceps

Doughanline's ziurki crantys aligatorianus
ngsru enyplis
Asidolder / vojamos



- Strong grasping forceps - designed for ERCP stent removal
- Rat tooth and alligator jaw for firm grip
- Rotatable for easy positioning




Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Features	Opening Width	Area Compatibility
FG-44NR-1	026790	1	2.8 mm	1800 mm	Rotatable	6.9 mm	1

5.1
Endoskopo kanalas
Darbinis ilgis
Rotuojama
Asidarymo plotis
stentų šalinimui

Product-Information-Sheet

Daughan Lin's iroclimo rinklins

Reusable 8.5 Fr. Conventional Stent Insertion Kit



- Radiopaque tip
- Excellent insertability
- Kit contains a guide catheter and a pusher tube



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Stent Size	Area Compatibility
MAJ-508	026517	1	3.2 mm	1600 mm	8.5 Fr.	①

Daughan Lin's iroclimo rinklins *Darbinis ilgis* *stentis dyelis*

Reusable 10 Fr. Conventional Stent Insertion Kit



- Radiopaque tip
- Excellent insertability
- Kit contains a guide catheter and a pusher tube




Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Stent Size	Area Compatibility
MAJ-255	026585	1	3.7 mm	1600 mm	10 Fr.	①

1
Darbinis ilgis *stentis dyelis*

Product-Information-Sheet

Danğhar tiri kanüle

Conventional Reusable Cannulae




Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Tip Shape	Distal Tip ø	Compatible Guidewire	brand_name	Area Compatibility
PR-304Q	N1078730	1	2.2 mm	1950 mm	Standard type	4 Fr.	0.035"	Star Tip, X-Press	①

Tin kama endoskopo kanale
Dan bini ilgis

Conventional Reusable Cannulae

Danğhar tiri kanüle



Article Name	Article number	Quantity	Min. Working Channel ø	Working Length	Tip Shape	Distal Tip ø	Compatible Guidewire	brand_name	Area Compatibility
PR-310Q	N1079130	1	2.2 mm	1950 mm	Long, tapered	3.5 Fr.	0.025"	Star Tip, X-Press	①

Endoskopo kanale
Dan bini ilgis
suderinanna su reddis

Product-Information-Sheet



Daughan's sarvas

Reusable Oval Electrosurgical Snares



- Wide-opening, oval snare
- SD-9U-1, -12U-1, and SD-310Z-25: Thick, braided wire (0.47 mm)
- SD-11U-1 and -13U-1: Conventional braided wire (0.43 mm)
- Outer plastic sheath to be combined with snare wire and SD handle (MH-264)


sudirinama sar ranlung



Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	For Use With	Area Compatibility
MAJ-15 8.3	SD-9U-1 to -13U-1	026008	1 sheath <i>1 sarvas</i>	2.8 mm	2300 mm	MAJ-16, -18, -20, or -22, MH-264	



Daughan's kitpos

Reusable Oval Electrosurgical Snares



- Wide-opening, oval snare
- SD-9U-1, -12U-1, and SD-310Z-25: Thick, braided wire (0.47 mm)
- SD-11U-1 and -13U-1: Conventional braided wire (0.43 mm)
- Outer plastic sheath to be combined with snare wire and SD handle (MH-264)


sudirinama sar ranlung



Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	Wire ø	Loop ø (A)	For Use With	Area Compatibility
MAJ-22 8.2	SD-13U-1	026022	1 wire	2.8 mm	2300 mm	0.43 mm	15 mm	MAJ-15 and MH-264	
MAJ-18 8.1	SD-11U-1	026017	1 wire	2.8 mm	2300 mm	0.43 mm	25 mm	MAJ-15 and MH-264	

Daughan's ranlung

Reusable Oval Electrosurgical Snares



- Wide-opening, oval snare
- SD-9U-1, -12U-1, and SD-310Z-25: Thick, braided wire (0.47 mm)
- SD-11U-1 and -13U-1: Conventional braided wire (0.43 mm)
- Outer plastic sheath to be combined with snare wire and SD handle (MH-264)



Autoblastarjams

Article Name	Article number	Quantity	For Use With
MH-264 8.4	027004	1 handle	SD-9U-1 to -13U-1

Product-Information-Sheet

Daugardini's sarvas

Single-Use Injection Needles with Reusable Metal Sheath



- Reusable outer metal-coil sheath used with single-use inner needles
- Ergonomic handle with positive ?click? action
- Metal-coil sheath enhances mucosal penetration
- Outer metal-coil sheath (MAJ-74, -67, -69) to be combined with inner needle (MAJ-75, -76, -68, -71, -72, -73, -70)






Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	For Use With	Area Compatibility
MAJ-67	NM-4 to 7L-1	026991	1 sheath	○ 2.8 mm	1650 mm	MAJ-68, -71, -72, -73	 

9.2 *Daugardini's sarvas* *Darbin's ilgi's*

Single-Use Injection Needles with Reusable Metal Sheath



- Reusable outer metal-coil sheath used with single-use inner needles
- Ergonomic handle with positive ?click? action
- Metal-coil sheath enhances mucosal penetration
- Outer metal-coil sheath (MAJ-74, -67, -69) to be combined with inner needle (MAJ-75, -76, -68, -71, -72, -73, -70)





Article Name	Component Of	Article number	Quantity	Min. Working Channel ø	Working Length	For Use With	Area Compatibility
MAJ-69	NM-4U-1	026475	1 sheath	○ 2.8 mm	2300 mm	MAJ-70	 

9.1 *Daugardini's sarvas* *Darbin's ilgi's*

Product-Information-Sheet

Daughantini avarinio lito hipobrius vankana

Reusable Emergency Lithotripter



- To assist in situations when baskets become impacted
- Set consists of a metal-coil sheath (MAJ-403) and a handle



Article Name	Article number	Quantity	Working Length	Features	Area Compatibility
BML-110A-1	026574	1 set	900 mm	Handle and sheath	①

10.1 *Daughantini avarinio lito hipobrius sarvas*

Reusable Emergency Lithotripter



- To assist in situations when baskets become impacted
- Set consists of a metal-coil sheath (MAJ-403) and a handle



Article Name	Article number	Quantity	Working Length	Features	Area Compatibility
MAJ-403	026575	1 sheath	900 mm	Sheath	①

10.2 *Larbinis ilgis*

Color Codes for Minimum Channel Diameters

The following color codes provide orientation to match EndoTherapy devices with the right Olympus channel diameter.

Minimum Channel Diameter	Color Code
1.2 mm/1.5 mm	○ White
1.7 mm	● Violet
2.0 mm/2.2 mm	● Blue
2.6 mm	● Green
2.8 mm/3.2 mm	● Yellow
3.7 mm/4.2 mm/6.0 mm	● Orange
5.0 mm/5.5 mm	● Pink

Symbols for Device Features

	Single use only		This device is sterile packed
	This device is designed for Refined Reprocessing (reusable, autoclavable)		This device is autoclavable
	This device has a luer-lock connector or cleaning adaptor for flushing		This device allows for the instillation of contrast medium
	This device is visible under fluoroscopic control		This device can be used in combination with electrosurgical devices
	This device is guidewire compatible		This device is compatible with a 0.018" guidewire
	This device is compatible with a 0.021" guidewire		This device is compatible with a 0.025" guidewire
	This device is compatible with a 0.035" guidewire		This device is compatible with a 0.038" guidewire

Symbols for Medical Fields

	For bronchoscopy		For upper GI endoscopy
	For biliary and pancreatic procedures		For lower GI endoscopy

OLYMPUS

Your Vision, Our Future

CONTENTS

01 Information & Service

Olympus – Your Vision, Our Future
 Olympus Medical Systems
 Social and Environmental Responsibility. . .
 Olympus EndoTherapy
 Guaranteed Compatibility between Original D
 and Ancillaries
 Quality Assurance
 Direct Delivery for Improved Customer Satisfac
 Olympus – Your Reliable Service Partner . .
 Olympus Nomenclature.
 Endoscope Compatibility – Color Coding . .

02 Biopsy Forceps

Single-Use Devices

Biopsy Forceps
 Hot Biopsy Forceps

Reusable Devices

SwingJaw Biopsy Forceps
 Biopsy Forceps
 Hot Biopsy Forceps

03 Polypectomy, EMR & ESD

Single-Use Devices

Electrosurgical Snares
 Hot Biopsy Forceps
 Grasping Forceps
 Ligation Devices
 Spray Catheter
 Distal Attachments