

## SonoSurg GE Transducers

N1806330 Transducer "SonoSurg T2L-GE", with connecting cable, with HF connector, for ultrasonic aspiration in open and endoscopic surgery



N1806430 Transducer "SonoSurg T2L-GE-C", with HF connector, for ultrasonic aspiration in open and endoscopic surgery, requires detachable connecting cable MAJ-1121



N1021500 Cable "MAJ-1121", for transducer T2L-GE-C




## SonoSurg GE Sheaths and Probes

N1807230 SonoSurg GE sheath "T3330", tapered, for endoscopic surgery



N1806830 Probe "MH-125", for sheath T3330



N1807130 SonoSurg GE sheath "T3320", 15° angled, for open surgery



N1806730 Probe "MH-113", for sheath T3320



N1807030 SonoSurg GE sheath "T3310", straight, for open surgery



N1806630 Probe "MH-112", for sheath T3310




N1806930 Attachment "T3300", for SonoSurg GE sheaths, 3 pcs.



## Accessories

N1807930 Tubing "MAJ-1099", set, for suction/irrigation, 6 pcs.



*Plavimol/sarbinol  
zarnebis*

## System Chart

